

Jerudong Park Medical Centre



**Organization Accredited Joint
Commission International**



LABORATORY HANDBOOK

DOCUMENT NUMBER:	LAB-HB-001/Version 7
EFFECTIVE DATE:	20/08/2024
AUTHORED BY:	Hazimah Hanis Hamzah (Scientific Officer, Lab Quality Coordinator) Lim Kai Wern (Scientific Officer) Rebecca Thien (Scientific Officer) Racheal Tang Lee Chau (Scientific Officer) Pg. Muhd. Nazri Pg. Mohd. Juanda (Scientific Officer) Linaa Farhana Abdillah (Lab Technician) Afiqkhazmiezul Ramdani (Lab Technician)
APPROVED BY:	Dr Hj Mohamad Hj Kassim (Consultant Microbiologist-cum-Laboratory Manager)

REVISION HISTORY		
Version #	Date Revised	Description
1	01/03/2016	Handbook created
2	14/10/2019	Added Urgent Request, Amended Report and Outsourced Lab Tests policy. Added blood culture collection, Added list of urine preservatives, Added new statements in Transportation of Specimen, Updated Critical Value List, Updated STAT tests with TAT, Removed Andrology and Embryology section, Rearranged test catalogue in alphabetical sauce. <u>Added new tests:</u> Vit D, Albumin (Fluid/CSF), Alcohol (Blood/Urine), Amylase (Fluid/Urine), Blood Gas, ACR, CCR, Glucose (Urine/Fluid), Cortisol Urine, Phenytoin, Phenobarbital, UIBC, TIBC, PIVKA, Theophylline, Trop-I, Trop-I (Risk), Triglycerides Fluid, Phosphate (Urine), Acinetobacter screening, Respiratory Viral Panel (CR), Respiratory Virus screening, A-HBc IgM, A-HBe, HBeAg, Hep B confirmatory, Hep C confirmatory, HIV Ab Confirmatory <u>Removed tests:</u> IgA, IgG CSF, IgM, Bleeding Time, Semen C/S, Chlamydia trachomatis Antigen, DNP, Respiratory Adenovirus, RSV
3	05/07/2021	Updated Organisation Chart, Added Add-On Request policy, Added new statements under Specimen Collection for Laboratory Testing, Added new statements under Transportation of Specimen, Updated critical result list <u>Added new tests:</u> Cortisol (Urine), Cortisol (24hr urine), Cryoglobulin test, SARS-CoV-2 IgM, SARS-CoV-2 IgG, Fungal Microscopy and Culture, Salmonella typhi IgM/IgG/ <u>Removed tests:</u> PIVKA, Glucose (Urine), Infectious Mononucleosis/ <u>Revised test info:</u> Cortisol, Urine Drug Screening, HbA1c, Microalbumin, RSCS, Faecal Occult Blood, C.difficile toxin, Faecal Adenovirus Antigen, Faecal Rotavirus Antigen, ZN stain
4	24/08/2022	Added safe disposal under specimen collection, Updated 24 hour urine preservatives, Added nasal/nasopharyngeal specimen collection, Updated diagram for specimen labelling, Updated critical result list, Updated list of STAT tests, Added Test Catalogue (Send-out), Moved Test Catalogues and List of STAT tests as annex documents

		<p><u>Added new tests:</u> CT/NG, CSF Lactate, CSF LDH, Cyclosporine, Ethanol, Ethanol (Urine), Fluid Albumin, Fluid Amylase, Fluid Lactate, Fluid LDH, Fluid Sodium, Fluid Triglycerides, Gastrointestinal PCR, Meningitis/Encephalitis PCR, Tacrolimus</p> <p><u>Removed tests:</u> Benzodiazepines, High Vaginal Swab C/S, TIBC, UIBC</p> <p><u>Revised test info:</u> Widal/Weil Felix, ESR, Direct Coombs Test, Cryoglobulin, Blood Culture, Antibody Screening, Rubella IgG, Rejection Criteria, Method on majority of Biochemistry tests</p>
5	27/01/2023	<p>Updated specimen requirement for urine drug to 2 urine specimen (page 14). Updated urine collection procedure (using vacuum). Update blood culture collection procedure.</p> <p><u>Removed tests:</u> COVID-19 IgM</p> <p><u>Revised test info:</u> Triglycerides, HDL Cholesterol, Cholesterol, Lipid Profile, Urine Drug Screen</p>
6	20/08/2023	<p>Updated and included copies of laboratory request forms. Removed ESR tube from list of specimen. Removed uric acid from list of urine that requires preservative</p> <p>Added “Leukodepleted packed cells” under blood products page 24.</p> <p><u>Added new tests:</u> NT-ProBNP, Respiratory Panel, C3, C4, Anti-CCP, Anti-Thyroglobulin, Anti-Thyroid Peroxidase, Thyroid Receptor Antibody, Thyroglobulin, Insulin, Parathyroid Hormone, Urine Glucose, Urine Ketone, HBV DNA Viral Load, Urea (Urine 24 hours)</p> <p><u>Removed tests:</u> COVID-19 IgG, BNP</p> <p><u>Added new outsourced tests:</u> NICE Basic, NICE Extended, Parathyroid Hormone</p> <p><u>Revised test info:</u> Widal/Weil Felix, Typhoid antibody, Rickettsial serology, Uric Acid, Amylase, Cryoglobulin Test, outsourced tests</p>
7	17/09/2024	<p>Added Customer Feedback Process. Updated laboratory request form. Updated Sample Transportation. Updated Blood Product request codes.</p> <p>Updated Transfusion Reaction sample requirements. Remove requirement for Blood Bank Request Form for emergency release case.</p> <p>Updated list of STAT tests. Change Blue Top volume requirement to “Up to indicated line on tube”.</p> <p><u>Added new tests:</u> Thyroid Autoantibodies, Fluid Chloride, Fluid Creatinine, Fluid Potassium, Fluid Urea, MTB/RIF</p> <p><u>Removed tests:</u> Glucose Tolerance Tests (3 Specimen), Immunoglobulin G, Thyroid Function Test 2 (TFT2)</p> <p><u>Revised test info:</u> Arterial/Venous Blood Gas, Diabetic Studies, Glomerular Filtration Rate. Preservative requirement for urine calcium and creatinine</p>

TABLE OF CONTENTS

INTRODUCTION	1
ORGANISATION CHART	2
LABORATORY OPENING HOURS.....	3
CUSTOMER FEEDBACK PROCESS.....	4
REQUESTS FOR LABORATORY TESTING	5
LABORATORY REQUEST FORM / ONLINE REQUEST	5
URGENT REQUEST	5
ADD-ON REQUEST	5
SPECIMEN COLLECTION FOR LABORATORY TESTING	9
BLOOD COLLECTION.....	9
BLOOD CULTURE COLLECTION	15
URINE COLLECTION	16
URINE DRUG SPECIMEN COLLECTION.....	19
STOOL COLLECTION	20
SPUTUM COLLECTION	20
NASOPHARYNGEAL SWAB COLLECTION	21
NASAL SWAB COLLECTION	21
SPECIMEN COLLECTION FOR FUNGUS	21
SURGICAL PATHOLOGY TISSUE COLLECTION	22
SUREPATH PAP TEST COLLECTION	22
SPECIMEN LABELLING	23
SPECIMEN HANDLING AND TRANSPORT	24
TRANSPORTATION OF SPECIMEN	24
REJECTION OF SPECIMEN	25
PRECIOUS SPECIMEN	25
LABORATORY REPORTING	26
AMENDED REPORT.....	26
OUTSOURCED LABORATORY TESTS	26
CRITICAL RESULTS	26
CRITICAL RESULT LIST	27
BLOOD TRANSFUSION POLICY	28

CROSSMATCH REQUEST..... 28

INVESTIGATION OF ADVERSE TRANSFUSION REACTION 29

EMERGENCY RELEASE OF BLOOD UNITS 30

ANNEX A. TEST CATALOGUE (IN-HOUSE).....

ANNEX B. LIST OF STAT TESTS.....

ANNEX C. TEST CATALOGUE (OUTSOURCE)

ANNEX D. TEST PACKAGES.....

JPMC PACKAGES.....

BNSRC PACKAGES.....

TBCC PACKAGES

GJPMC PACKAGES

INTRODUCTION

OVERVIEW OF JPMC LABORATORY

As part of the diagnostic services of Jerudong Park Medical Centre, the Department of Laboratory Services is committed in providing high quality laboratory services relevant to medicine and disease surveillance to physicians and health care providers for optimum patient care.

Laboratory testing is performed 24 hours a day, seven days a week. Certain tests are referred to approved and accredited reference laboratories.

The Department of Laboratory Services operates the following divisions:

- I. Biochemistry and Immunology
- II. Blood Bank and Transfusion
- III. Haematology
- IV. Clinical Microbiology
- V. Virology and Serology
- VI. Specimen Receiving and Management
- VII. Phlebotomy
- VIII. Blood Donation

MISSION STATEMENT

Our department mission is to provide high-quality, timely, reliable and accurate services, using innovative technology, to meet the needs and requirement of our customers.

VISION STATEMENT

Our department vision is to be the premier provider of medical laboratory services in our community.


QUALITY ASSURANCE

Our quality management system adheres to our customer requirement, ISO 15189 and JCI accreditation standards to ensure the accuracy and quality of our test results. Our laboratory also participates in external quality assessment (EQA) programs, such as College of American Pathologists (CAP) and Royal College of Pathologist of Australasia (RCPA).

Laboratory only claims conformity to ISO15189 for tests listed in the scope of accreditation.

Refer to changes highlighted yellow

ORGANISATION CHART



HOS – Head of Section
 MSO – Medical Scientific Officer
 MLT – Medical Lab Technologist
 LT- Lab Technician

Refer to changes highlighted yellow

LABORATORY OPENING HOURS

*Laboratory will be in operation 24/7 for in-patient and urgent outpatient requests

MONDAY TO THURSDAY	08:00AM - 12:00PM 01:00PM - 04:00PM
FRIDAY	08:00AM - 12:00PM 02:00PM - 05:00PM
SATURDAY	08:00AM - 12:00PM
SUNDAY & PUBLIC HOLIDAY	CLOSED

PHLEBOTOMY SERVICE (JPMC)

MONDAY TO THURSDAY	07:30AM - 12:00PM 01:00PM - 03:30PM
FRIDAY	07:30AM - 12:00PM 02:00PM - 04:30PM
SATURDAY	08:00AM - 12:00PM
SUNDAY & PUBLIC HOLIDAY	CLOSED

PHLEBOTOMY SERVICE (PJSC)

MONDAY TO THURSDAY	07:30AM - 12:00PM 01:00PM - 04:30PM
FRIDAY	07:30AM - 12:00PM 02:00PM - 04:30PM
SATURDAY, SUNDAY & PUBLIC HOLIDAY	CLOSED

Refer to changes highlighted yellow

For enquiries, we can be reached at:

Location	Level 1, Jerudong Park Medical Centre
Telephone	+6732611433 (see extension below)
Email	lab.reception@jpmc.com.bn

Section	Extension
Laboratory Manager	2201 / 2382
Administration	2415
Quality Assurance	2416
Phlebotomy (JPMC)	2516
Phlebotomy (PJSC)	(8) 1575
Central Specimen Receiving and Management	2127
Biochemistry and Immunology	2407
Haematology	2408
Virology and Serology	2322
Microbiology	2436
Blood Bank and Transfusion	2320
Blood Donation	2414
Laboratory Store	2381
Point-of-Care (POCT)	2407

CUSTOMER FEEDBACK PROCESS

For any feedback about our laboratory services, please contact us using the information provided above.

1. **Contact Us:** Reach out to us via email, phone, customer feedback form, IR, or in person.
2. **Feedback Assessment:** Your feedback will be reviewed and assessed by our team.
3. **Investigation:** If your feedback requires further investigation, we will take appropriate steps to address your concerns.
4. **Feedback Response:** We will provide a timely response to your feedback, including the results of any investigation if necessary.

REQUESTS FOR LABORATORY TESTING

LABORATORY REQUEST FORM / ONLINE REQUEST

- ✓ It is essential that under all circumstances the request (online and /or via hard copy) be filled out COMPLETELY.
- ✓ The laboratory shall not perform any examination on a specimen without a request form (online or hard copy). **Verbal request is not accepted as an alternative to a request form.**
- ✓ All request details must be filled in and the form MUST clearly show to whom the results are to be reported / communicated to.
- ✓ The ordering physician is responsible for ensuring the following information is entered LEGIBLY and CORRECTLY on the request form:
 - a. Patient's name
 - b. Patient's MRN/PRN
 - c. Patient's identity card no. (IC No)
 - d. Date and year of birth
 - e. Sex of patient
 - f. Ordering location
 - g. Date and time of specimen collection
 - h. Specimen type and/or source
 - i. Name and signature of ordering physician
 - j. Relevant clinical data and diagnosis
 - k. Test requests
- ✓ Manual order must be filled using the "LABORATORY REQUEST FORM" for general tests and "MICROBIOLOGY REQUEST FORM" for Microbiology tests.
- ✓ Blood product request form must be filled using the "BLOOD BANK REQUEST FORM".

URGENT REQUEST

- ✓ Urgent request is marked as appropriate in the request form and will be given priority over routine requests. Urgent request should only be used for results that are required as soon as possible for immediate patient management.
- ✓ Only Short Turn-Around-Time (STAT) tests are processed as Urgent. STAT tests are listed in **Annex B**.

ADD-ON REQUEST

- ✓ Ordering location must send a new requisition form indicating it is an add-on request. Laboratory accepts add-on request if specimen is acceptable for additional testing (i.e. stability, volume). If specimen is no longer suitable, request will be rejected and requesting location informed.

Refer to changes highlighted yellow



Laboratory Numbers

PRN No.: _____

Name: _____

IC/Passport No.: _____

Sex: _____ DOB: _____

Nationality: _____ Tel: _____

Please fill in or stick patient's sticker here

LAB REQUEST FORM

URGENT ROUTINE

REQUESTOR DETAILS

Requesting Doctor's Name: _____

Clinic / Ward: _____

Tel Number: _____

Appointment Date: _____
(AM/PM)

SAMPLE DETAILS

Sampling Date & Time: _____

24 hrs Urine Total Volume: _____ mL

Sample Taken By: _____

Fasting Non-Fasting

Patient Weight: _____ Kg

Patient Height: _____ Cm
(For creatinine clearance/Paediatric eGFR (<18 y.o.))

SAMPLE TYPE

Blood

Urine

Fluids: _____

Other: _____

RELEVANT CLINICAL INFORMATION

Drug Therapy: _____ Last Dose: _____
(Please note anticoagulant therapy when ordering coagulation test)

Other Relevant Clinical Information: _____ Date & Time: _____

TEST REQUESTS

PACKAGES	SINGLE TESTS	
<p>BNSRC PACKAGES</p> <p><input type="checkbox"/> NSRC1 <input type="checkbox"/> NSRC2 <input type="checkbox"/> NSRC3</p> <p>GJ PACKAGES</p> <p><input type="checkbox"/> ANGIOJ <input type="checkbox"/> GT5GJ <input type="checkbox"/> RFTGJ</p> <p><input type="checkbox"/> GJFRL <input type="checkbox"/> GT7GJ</p> <p><input type="checkbox"/> GTIGJ <input type="checkbox"/> OHS</p> <p>TBCC PACKAGES</p> <p><input type="checkbox"/> BRCP <input type="checkbox"/> LICP <input type="checkbox"/> NPC2</p> <p><input type="checkbox"/> COCP <input type="checkbox"/> LUCP <input type="checkbox"/> PACP</p> <p><input type="checkbox"/> GCCP <input type="checkbox"/> LYMCP <input type="checkbox"/> PRCP</p> <p><input type="checkbox"/> GYCP <input type="checkbox"/> NCCP <input type="checkbox"/> TPTP1</p> <p><input type="checkbox"/> HB6 <input type="checkbox"/> NPC1 <input type="checkbox"/> TPTP2</p> <p>PROFILES</p> <p><input type="checkbox"/> ELECTROLYTES (ELY)</p> <p><input type="checkbox"/> ELECTROLYTES, CREATININE AND UREA (ECU)</p> <p><input type="checkbox"/> CARDIAC ENZYMES (CAR)</p> <p><input type="checkbox"/> COAGULATION PROFILE (COP)</p> <p><input type="checkbox"/> LIVER FUNCTION TEST (LFT)</p> <p><input type="checkbox"/> RENAL FUNCTION TEST (RFT)</p> <p><input type="checkbox"/> THYROID BLOOD GAS (TFT)</p> <p><input type="checkbox"/> GLOMERULAR FILTRATION RATE (GFR)</p> <p><input type="checkbox"/> LIPID STUDIES (LIP)</p> <p><input type="checkbox"/> IRON STUDIES (IS)</p> <p><input type="checkbox"/> URINE DRUG SCREENING (UDRUG)</p> <p><input type="checkbox"/> ARTERIAL BLOOD GAS (GAS)</p> <p><input type="checkbox"/> VENOUS BLOOD GAS (VBG)</p> <p><input type="checkbox"/> PROTEIN/CREATININE RATIO (PCRTR)</p> <p><input type="checkbox"/> ALBUMIN/CREATININE RATIO (ACR)</p>	<p>BIOCHEMISTRY</p> <p><input type="checkbox"/> ALBUMIN</p> <p><input type="checkbox"/> ALP</p> <p><input type="checkbox"/> ALT</p> <p><input type="checkbox"/> AMIKACIN - PEAK</p> <p><input type="checkbox"/> AMIKACIN - TROUGH</p> <p><input type="checkbox"/> AMMONIA</p> <p><input type="checkbox"/> AMYLASE</p> <p><input type="checkbox"/> AST</p> <p><input type="checkbox"/> BICARBONATE</p> <p><input type="checkbox"/> BILIRUBIN DIRECT</p> <p><input type="checkbox"/> BILIRUBIN TOTAL</p> <p><input type="checkbox"/> C3</p> <p><input type="checkbox"/> C4</p> <p><input type="checkbox"/> CALCIUM</p> <p><input type="checkbox"/> CARBAMAZEPINE</p> <p><input type="checkbox"/> CHLORIDE</p> <p><input type="checkbox"/> CHOLESTEROL</p> <p><input type="checkbox"/> CK</p> <p><input type="checkbox"/> CREATININE</p> <p><input type="checkbox"/> CRP</p> <p><input type="checkbox"/> DIGOXIN</p> <p><input type="checkbox"/> ETHANOL</p> <p><input type="checkbox"/> GENTAMICIN - PEAK</p> <p><input type="checkbox"/> GENTAMICIN - TROUGH</p> <p><input type="checkbox"/> GGT</p> <p><input type="checkbox"/> GLUCOSE - FASTING</p> <p><input type="checkbox"/> GLUCOSE - RANDOM</p> <p><input type="checkbox"/> GLYCATED HAEMOGLOBIN (HBAIC)</p> <p><input type="checkbox"/> IRON</p> <p><input type="checkbox"/> LACTATE</p> <p><input type="checkbox"/> LDH</p> <p><input type="checkbox"/> MAGNESIUM</p> <p><input type="checkbox"/> MICROALBUMIN</p> <p><input type="checkbox"/> OSMOLALITY - BLOOD</p> <p><input type="checkbox"/> OSMOLALITY - URINE</p> <p><input type="checkbox"/> PHENOBARBITAL</p>	<p><input type="checkbox"/> PHENYTOIN</p> <p><input type="checkbox"/> PHOSPHATE</p> <p><input type="checkbox"/> POTASSIUM</p> <p><input type="checkbox"/> RHEUMATOID FACTOR</p> <p><input type="checkbox"/> SODIUM</p> <p><input type="checkbox"/> THEOPHYLLINE</p> <p><input type="checkbox"/> TOTAL PROTEIN</p> <p><input type="checkbox"/> TRANSFERRIN</p> <p><input type="checkbox"/> TRIGLYCERIDES</p> <p><input type="checkbox"/> UREA</p> <p><input type="checkbox"/> URIC ACID</p> <p><input type="checkbox"/> VALPROIC ACID</p> <p><input type="checkbox"/> VANCOMYCIN - PEAK</p> <p><input type="checkbox"/> VANCOMYCIN - TROUGH</p> <p>IMMUNOLOGY</p> <p><input type="checkbox"/> AFP</p> <p><input type="checkbox"/> ANTI-CCP</p> <p><input type="checkbox"/> ANTI-THYROGLOBULIN</p> <p><input type="checkbox"/> ANTI-THYROID PEROXIDASE</p> <p><input type="checkbox"/> B12</p> <p><input type="checkbox"/> BETA-HCG</p> <p><input type="checkbox"/> CA 125</p> <p><input type="checkbox"/> CA 153</p> <p><input type="checkbox"/> CA 19-9</p> <p><input type="checkbox"/> CEA</p> <p><input type="checkbox"/> CKMB</p> <p><input type="checkbox"/> CORTISOL</p> <p><input type="checkbox"/> DHEAS</p> <p><input type="checkbox"/> E2</p> <p><input type="checkbox"/> FERRITIN</p> <p><input type="checkbox"/> FOLATE</p> <p><input type="checkbox"/> FSH</p> <p><input type="checkbox"/> FT3</p> <p><input type="checkbox"/> FT4</p> <p><input type="checkbox"/> INSULIN</p> <p><input type="checkbox"/> LH</p> <p><input type="checkbox"/> METHOTREXATE</p>
	<p><input type="checkbox"/> NT-PROBNP</p> <p><input type="checkbox"/> PARATHYROID HORMONE</p> <p><input type="checkbox"/> PROCALCITONIN</p> <p><input type="checkbox"/> PROGESTERONE</p> <p><input type="checkbox"/> PROLACTIN</p> <p><input type="checkbox"/> TOTAL PSA</p> <p><input type="checkbox"/> TOTAL TESTOSTERONE</p> <p><input type="checkbox"/> THYROGLOBULIN</p> <p><input type="checkbox"/> TROPONIN I</p> <p><input type="checkbox"/> TSH</p> <p><input type="checkbox"/> TSH RECEPTOR ANTIBODY</p> <p><input type="checkbox"/> 25-OH VITAMIN D TOTAL</p> <p>HAEMATOLOGY</p> <p><input type="checkbox"/> ACTIVATED PARTIAL THROMBOPLASTIN TIME</p> <p><input type="checkbox"/> APTT 50% CORRECTION</p> <p><input type="checkbox"/> APTT RATIO (APTT INCLUDED)</p> <p><input type="checkbox"/> BLOOD FILM</p> <p><input type="checkbox"/> D-DIMER</p> <p><input type="checkbox"/> ESR</p> <p><input type="checkbox"/> FBC</p> <p><input type="checkbox"/> FIBRINOGEN</p> <p><input type="checkbox"/> HAEMOGLOBIN</p> <p><input type="checkbox"/> MALARIA SCREEN</p> <p><input type="checkbox"/> MICROFILARIA SCREEN</p>	<p><input type="checkbox"/> PROTHROMBIN TIME</p> <p><input type="checkbox"/> PROTHROMBIN 50% CORRECTION</p> <p><input type="checkbox"/> RETICULOCYTE</p> <p>VIROLOGY/SEROLOGY</p> <p><input type="checkbox"/> DENGUE SEROLOGY</p> <p><input type="checkbox"/> FLU/RSV</p> <p><input type="checkbox"/> HEP A AB</p> <p><input type="checkbox"/> HEP B CORE AB</p> <p><input type="checkbox"/> HEP B SURFACE AB</p> <p><input type="checkbox"/> HEP B SURFACE AG</p> <p><input type="checkbox"/> HEP C AB</p> <p><input type="checkbox"/> HIV SCREENING</p> <p><input type="checkbox"/> RPR</p> <p><input type="checkbox"/> SYPHILIS AB</p> <p>BLOOD TRANSFUSION</p> <p><input type="checkbox"/> ANTIBODY SCREENING</p> <p><input type="checkbox"/> BLOOD GROUP</p> <p><input type="checkbox"/> BLOOD GROUP AND HOLD</p> <p><input type="checkbox"/> DIRECT COOMBS TEST</p>

ADDITIONAL TESTS REQUESTED

Doctor's Signature: _____ Date: _____

FOR LAB USE			
BLOOD	PLAIN	EDTA	GREY
	GREEN	BLUE	
URINE	Received By and Date:		
SWAB			
FLUID			
OTHERS			

Refer to changes highlighted yellow



MICROBIOLOGY REQUEST FORM

PRN No.: _____
 Name: _____
 IC/Passport No.: _____
 Sex: _____ DOB: _____
 Nationality: _____ Tel: _____
Please fill in or stick patient's sticker here

Laboratory Numbers

THIS SECTION MUST BE COMPLETED TO AVOID DELAYS IN SAMPLE PROCESSING

Priority: (Routine unless otherwise stated) <input type="checkbox"/> URGENT All tests: <input type="checkbox"/> Yes <input type="checkbox"/> No If no, specify tests: _____ <input type="checkbox"/> Scheduled appointment date: _____ (Date) (Time)	SAMPLE COLLECTION Date of collection: _____ Time of collection: _____ Taken by: _____ Checked by: _____	ORDERING PHYSICIAN Name: _____ Ward/Clinic: _____ Signature of Physician: _____ Date Signed: _____
---	--	---

DIAGNOSIS	CURRENT MEDICATION/ DATE & TIME OF LAST DOSE

URINE <input type="checkbox"/> Routine Urinalysis <input type="checkbox"/> Urinalysis + Culture For culture, list current antibiotics: _____	STOOL <input type="checkbox"/> Microscopy (ova, cyst, parasites) <input type="checkbox"/> Microscopy + Culture <input type="checkbox"/> Faecal Occult Blood / Faecal Immunoassay Test <input type="checkbox"/> C.difficile Test	BLOOD CULTURE <input type="checkbox"/> Aerobic / Anaerobic Sites: _____ <input type="checkbox"/> Paediatric <input type="checkbox"/> Myco / Lytic	CSF/FLUID Source: _____ <input type="checkbox"/> Cell count with differential <input type="checkbox"/> Gram Stain <input type="checkbox"/> Culture <input type="checkbox"/> Cryptococcus (India Ink) <input type="checkbox"/> Others: _____
--	--	--	--

M/E (GRAM STAIN - IF APPLICABLE) & CULTURE Specimen Type: <input type="checkbox"/> Respiratory (Throat, Sputum, Bronchoalveolar lavage, Endotracheal aspirate, Nasopharyngeal aspirate) <input type="checkbox"/> Wound / Pus / Abscess Sites: 1. _____ 2. _____ <input type="checkbox"/> Eye <input type="checkbox"/> Ear <input type="checkbox"/> Nasal <input type="checkbox"/> Semen <input type="checkbox"/> Tip: _____ <input type="checkbox"/> Genital - Direct smear (Trichomonas & yeasts) and/or Gram stain & Culture: <input type="checkbox"/> High vaginal swab <input type="checkbox"/> Urethral (Male/Female) <input type="checkbox"/> Low vaginal swab <input type="checkbox"/> Penile <input type="checkbox"/> Group B Streptococcus (Antenatal) <input type="checkbox"/> MRSA Screening Sites: _____ <input type="checkbox"/> Acinetobacter Screening Sites: _____ <input type="checkbox"/> Sterility Test Sites: _____	MYCOLOGY <input type="checkbox"/> KOH prep (Direct Smear) <input type="checkbox"/> Culture Specimen Type: <input type="checkbox"/> Skin <input type="checkbox"/> Nail <input type="checkbox"/> Hair <input type="checkbox"/> Other site: _____	OTHER TESTS <input type="checkbox"/> AFB Smear & Culture Sites: _____ <input type="checkbox"/> Fungal Culture Sites: _____ <input type="checkbox"/> Widal <input type="checkbox"/> Weil Felix <input type="checkbox"/> Salmonella typhi IgG / IgM <input type="checkbox"/> Meningitis / Encephalitis Panel <input type="checkbox"/> Gastrointestinal Panel
---	--	--

ADDITIONAL NOTES / TESTS:

LABORATORY USE ONLY

Blood		Others	
Urine		Sputum	
Stool		BC	
Swab		Fluid	

Received by/time/date: _____

Refer to changes highlighted yellow



LABORATORY DEPARTMENT

LAB NO: _____

BLOOD BANK REQUEST FORM

Doctor's Name _____	MRN _____
Clinic / Ward _____	Name _____
Date _____	ID No. _____
Time _____	Sex _____
Signature _____	D.O.B _____

<p>SPECIMEN: 1 EDTA Blood Tube</p> <p>Coll'd by: _____ Rec'd by: _____ Date: _____ Date: _____ Time: _____ Time: _____</p>	<p>Clinical History:</p> <p>Hb (g/L): _____ PLT (10³/μL): _____</p>
---	--

BLOOD TRANSFUSION REQUEST

Please Check the , and indicate No. of Units or amount (mL)

Further Information

<input type="checkbox"/> Packed Cells _____	Blood Group (if known) _____	Rh(D) _____
<input type="checkbox"/> Whole Blood _____	For Surgery / Transfusion on _____	
<input type="checkbox"/> Fresh Frozen Plasma (FFP) _____		
<input type="checkbox"/> Cryoprecipitate _____		
<input type="checkbox"/> Random Platelet Units _____		
<input type="checkbox"/> Apheresis Units _____		

BLOOD TRANSFUSION RECORD

	Unit No.	Group/Rh	Expiry	Vol (mL)	Issued by	Taken by	Location	Date & Time	Comment
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									

THAWED FFP/CRYO
 Date and time thawed: _____
 FFP unit/s must be transfused within 24 hours after thawing.
 CRYO unit/s must be transfused within 6 hours after thawing.

Technologist Signature _____
 Date _____
 Time _____

Remarks: _____

NOTE:
 • Return blood to Blood Bank within 30 minutes if not used. • Blood will be held only until 8 am on the 5th calendar day after surgery or reserve date, unless Blood Bank is notified.

BLOOD TRANSFUSION LAB : Ext 2320

JPMC/BB/BBRF/VERO3/11.02.20

SPECIMEN COLLECTION FOR LABORATORY TESTING

Perform hand hygiene and wear PPE before collecting specimen from patients


- ✓ Laboratory test results are dependent on the quality of the specimen submitted. If there is any doubt or question regarding the type of specimen that should be collected, please call Laboratory for clarification.
- ✓ The person collecting the sample is responsible in identifying the patient (using at least TWO IDENTIFIERS: Patient Full Name and Date of Birth) prior to collection.
- ✓ Ensure patient is prepared as per test requirement (e.g. fasting status, specimen collected at timed intervals, medication status). Refer to **Annex A**. Test Catalogue for individual test requirement.
- ✓ When there is deviation to the documented collection procedure (e.g. patient not fasting for fasting glucose test, or blood specimen taken from hematoma site), the deviation must be recorded in the request form or communicated to the laboratory personnel.
- ✓ Ensure specimen collection and transportation supplies are not expired.
- ✓ After specimen collection, dispose all disposable PPE and consumables used into the yellow biohazard waste bag. Dispose sharp items (such as needles) into leak-proof and puncture resistant sharps waste container.






BLOOD COLLECTION

- ✓ Collect the blood specimen from a vein and avoid prolonged application of tourniquet (must be less than or within 1 minute) to prevent stasis and hemolysis.
- ✓ To Avoid Hemolysis:
 - Ensure alcohol has dried before venipuncture
 - Use the correct needle gauge when collecting blood sample for pediatric patients and for patient with difficult veins
 - Do not mix the blood tubes vigorously
 - Do not collect specimens from veins where administration of fluids will cause abnormal levels of electrolytes, glucose and drugs.
 - Avoid contamination from heparin locks for coagulation test
 - Collect the specimen into the proper tube or container using the correct sequence of draw
- ✓ The following order of draw is recommended when drawing multiple specimen during a single venipuncture:

Refer to changes highlighted yellow

PAEDIATRIC PATIENTS



CONTAINER	BLOOD CULTURE	LIGHT BLUE TOP	PURPLE MICROTAINER	GREEN MICROTAINER	YELLOW MICROTAINER
					
ADDITIVE	Bacterial Growth Medium	Sodium citrate	K ₂ EDTA	Lithium Heparin	Clot Activator
GENTLE MIX BY INVERSION	8-10 times	6-8 times	10 times	10 times	5 times
MINIMUM DRAW VOLUME	3 mL	Up to line	500µL	600µL	600µL

For GTP package, 1 x Lavender Top and 2 x Green Top preferred.


Refer to changes highlighted yellow

ADULT PATIENTS

CONTAINER	ADDITIVES	GENTLE MIX BY INVERSION*	MINIMUM DRAW VOL.
 <ul style="list-style-type: none"> BLOOD CULTURE (AEROBIC & ANAEROBIC BOTTLES) 	BACTERIAL GROWTH MEDIUM	8-10 times	8-10 mL
 <ul style="list-style-type: none"> LIGHT BLUE TOP 	SODIUM CITRATE	3-4 times	Up to line
 <ul style="list-style-type: none"> SST PLAIN/YELLOW TOP 	CLOT ACTIVATOR AND GEL FOR SERUM SEPARATOR	5 times	5 mL
 <ul style="list-style-type: none"> GREEN TOP 	LITHIUM HEPARIN	8-10 times	4 mL
 <ul style="list-style-type: none"> LAVENDER/PURPLE TOP 	K2 EDTA	8-10 times	4 mL
 <ul style="list-style-type: none"> ROYAL BLUE TOP 	K2 EDTA	8-10 times	6 mL
 <ul style="list-style-type: none"> GREY TOP 	FLOURIDE OXALATE	8-10 times	2 mL

ORDER OF DRAW

IMPORTANT Invert tubes immediately after collection by gentle inversion (see figure for example of 1 inversion). Improper mixing may cause clot in anti-coagulated blood. Do not shake tube as this may cause haemolysis.



Refer to changes highlighted yellow



Step 1 Apply the tourniquet about 4-6 inches above the venipuncture site. Do not leave tourniquet on for > 1 minute.



Step 2 The vein should be visible while applying the tourniquet. Locating the vein will help in determining the correct size of needle.



Step 3 Clean venipuncture site by 70% Alcohol swab in circular motion starting from the centre moving outward.



Step 4 Wait for alcohol to dry. Avoid touching cleaned site.



Step 5 Insert the needle gently.



Step 6 Once the needle hit the vein, the tube should be insert forward into the holder to puncture the stopper to collect the blood.



Step 7 When the blood begins to flow into the tube tourniquet should be release without moving the needle.



Step 8 The tube should filled until the vacuum is exhausted.
Note: it is important that the evacuated tube be filled completely because some additive tube are provided based on their full collection.



Step 9 Withdraw needle gently.



Step 10 Apply pressure on pricked site using sterile gauze.

Refer to changes highlighted yellow

Rejection Criteria	Main Causes	Corrective Action
Haemolysis	Wrong size of needle (Needle too small)	Use correct size of needle
	Venipuncture performed before alcohol is allowed to dry	Allow alcohol to dry thoroughly
	Prolonged tourniquet time	Release the tourniquet as soon as blood flow is established in the first tube. Limit tourniquet time to 1 minute or less
	Drawing blood on haematoma site	Select a different site. If different site not available, collect distally to the haematoma
	Difficult venipuncture/vein trauma/probing	The needle should be parallel to the vein. Enter at a 30° angle or less. Avoid probing
	Needle occlusion	Needle bevel may be positioned against the vein wall. Pull back slightly on the needle. Avoid rotating or changing the angle of the needle
	Drawing plunger back too forcefully in the syringe (syringe method)	Draw plunger back gently
	Air leakage around the needle or loss of vacuum in the tube	Make sure the needle is fitted securely on the syringe to avoid frothing
	Removing the needle from the vein with the tube still engaged (esp vacutainer method)	Remove tube first before removing the needle from vein
	Pushing on a syringe plunger too hard when expelling blood into a collection device	If there is loss of vacuum in the tube, safely remove the needle and gently transfer the syringe's contents by drizzling the blood down the side of the tube
	Underfilling tube (causes excessive blood-to-additive ratio)	Fill tubes to the correct volume
	Vigorously mixing the tube and frothing of sample	Use gentle inversion only
	Prolonged contact serum/plasma with cells – delayed centrifugation	Send sample to lab for centrifugation within 2 hours after collection
Exposure to high or low temperature before centrifugation	Do not place the tubes on a counter exposed to extreme temperature variation (e.g. next to autoclave, or directly touching ice pack – except for Ammonia sample). Never refrigerate the specimen before it is centrifuged	

Refer to changes highlighted yellow


	Centrifuging specimen before 30 minutes of collection (before clot formation for SST yellow tubes)	Allow the specimen tube to remain vertical in a rack for a minimum of 30 minutes
Haemolysis/Diluted specimen	Drawing blood from vascular/peripheral IV line	Pause the infusion for at least 2 minutes, flush the line with an adequate volume of saline, and discard enough waste blood to fully clear the flush prior to collecting lab blood specimen
Clotted	Difficult venipuncture/vein trauma/probing	The needle should be parallel to the vein. Enter at a 30° angle or less. Avoid probing
	Specimens not mixed well (poor inversion of the tube)	Mix well and gently invert specimen as soon as possible after collection to avoid clot formation

Refer to changes highlighted yellow

BLOOD CULTURE COLLECTION

IMPORTANT Must disinfect venipuncture site AND rubber top of blood culture bottles

- If using winged blood collection set, or using needle and syringe, or volume collected is less than recommended volume, inoculate aerobic blood culture bottle first if 1 set (one aerobic and one anaerobic) blood culture is required.
 - If 2 or more sets of blood culture is required, take each set from different sites. Indicate the site on the blood culture bottles and on the request form.
 - Collect specimen for culture ideally before administration of antimicrobial therapy.
-
- ◇ Choose a venipuncture site on the opposite extremity of an infusion.
 - ◇ Clean the venipuncture site with either 70% alcohol antiseptic pad or 2-4% chlorhexidine for at least 30 seconds.
 - ◇ Allow antiseptic or chlorhexidine to air dry completely. Do NOT palpate the site again.
 - ◇ Draw 16-20 mL of blood from adult patients, or 1-3 mL from paediatric patients.
 - ◇ Remove the flap covering the blood culture bottle and disinfect the rubber top with 70% alcohol pad prior to dispensing blood into them. Allow it to air dry.
 - ◇ Divide the blood by injecting the appropriate amount into each culture bottle, i.e. 8-10 mL of blood first into the aerobic bottle and then 8-10 mL into anaerobic bottle for adult patients.
 - ◇ Mix the blood adequately with the medium by inverting the bottles 8 to 10 times gently.
 - ◇ Discard syringes and needles in sharp container.
 - ◇ Label the bottles with patient identification details, the date and time of specimen collection. Do NOT put sticker over the bar code and lot number.




Do not put sticker or cover the bar code and lot number.

- ◇ Send specimen to the laboratory within 2 hours of collection. Specimen should be held at room temperature only. Do NOT refrigerate or freeze.

Refer to changes highlighted yellow

URINE COLLECTION

IMPORTANT Urine specimen for urinalysis and/or culture: 1) Date and time urine sample collected MUST be noted on the urine container, 2) Urine specimen stored at room temperature must be sent within 2-4 hours to the lab. Delayed delivery will result in sample rejection.



After urine specimen is collected in the urine container (Fig. 1), transfer urine to urine tubes (Fig.2 – prepare separate tubes for Biochemistry test and Microbiology test) according to the test requested.

Transfer as per instructions below:



Gently shake the sample

Partially raise the protective label (Do NOT remove it completely)

Insert the sample tube and make sure the needle penetrates the stopper of the urine tube. Keep the tube connected until it is full (end of flow)

Once tubes filled, remove the sample tube and fully re-stick the protective label

Shake the sample tubes 8-10 times, label the tubes with patient stickers and send to the Lab

- ◇ Do not send the urine container (Fig 1) to the laboratory: 1) unless urine volume is not sufficient to flow through the needle, or 2) for **CT/NG test**

Refer to changes highlighted yellow

Preservatives for random & 24 hours urine collection:

TEST	RANDOM URINE	24 HOUR URINE
Albumin/Creatinine Ratio	No preservative needed	20 to 30 mL of 6M HCL
Calcium	No preservative needed	20 to 30 mL of 6M HCL
Calcium/Creatinine Ratio	No preservative needed	20 to 30 mL of 6M HCL
Magnesium	No preservative needed	20 to 30 mL of 6M HCL
Microalbumin	No preservative needed	20 to 30 mL of 6M HCL
Phosphorus	No preservative needed	20 to 30 mL of 6M HCL
5 HIAA	NA	10 mL of 6M HCL (Send Out)
Metanephrines	NA	10 mL of 6M HCL (Send Out)
Vanillyl Mandelic Acid Homovanillic Acid	NA	20 mL of 6M HCL (Send Out)

- A) **Routine or random Specimen:** The patient is given a collection container and instructed to collect a midstream specimen in the container.
- B) **First voided Specimen:** First morning specimen or 8-hour specimen. The patient should be instructed to collect the specimen immediately after a night's sleep. The bladder is preferably emptied before sleep.
- C) **2-Hours postprandial Specimen:** The patient should be instructed to void shortly before consuming concentrated solution of 50g or 75g glucose or routine meal, and to collect a specimen 2 hours after eating.

Refer to changes highlighted yellow

D) **Clean-catch midstream Specimen:**

Patient Instructions: Male

- ◇ Wash hands thoroughly with soap and water and dry with a paper towel.
- ◇ Pass the initial portion of urine into the toilet bowl. Collect a portion of the remaining urine into a sterile collection container.

Patient Instructions: Female

- ◇ Wash hands thoroughly with soap and water and dry with a paper towel.
- ◇ Spread the labia and keep them apart.
- ◇ Cleanse the urethral meatus from front to back.
- ◇ Pass the initial portion of urine into the toilet bowl. Collect a portion of the remaining urine into a sterile collection container. Avoid contact of container with the legs, vulva, or clothing.

Transport the specimen to the laboratory immediately or refrigerate if transport is delayed.

E) **24-hour timed urine collection**

- ◇ Depending on the test, bring 24-hour container to laboratory prior to urine collection for it to be filled with respective preservative.
- ◇ To complete a 24 hour period, start and end collection at approximately the same time in the morning.
- ◇ Do not void directly into container.
- ◇ Collect each specimen in a disposable clean plastic or paper cup and carefully pour into the 24-hour container to avoid splatter or spillage.
- ◇ Patient should be instructed to avoid fecal contamination of the specimen.
- ◇ Refrigerate specimen during and after collection.

Day 1: Discard the first morning specimen (make sure you completely empty the bladder). Record on label: START DATE AND TIME. Begin collecting **ALL** subsequent specimen for the next 24 hours.

Day 2: Collect the last urine specimen exactly 24 hours after the start time, and then **STOP** Collection. Record on label: STOP (FINISH) DATE AND TIME.

- ◇ Tighten lid securely. Keep upright. Transport in a refrigerated container with requisition form as soon as possible after completion.
- ◇ If the amount of specimen exceeds the 24-hour urine container, collect the remaining specimen in a sterile collection container and label as "24-HOUR URINE #2". Note on the request form that 2 containers submitted.

Refer to changes highlighted yellow

- ◇ Dietary restrictions are required before and during the collection period for some tests (For example: Vanillyl Mandelic Acid test). Normal fluid intake is allowed during the collection period.

For Malay version of patient instruction, please call extension 2127

URINE DRUG SPECIMEN COLLECTION



Picture 1 Unused urine container



Picture 2 Properly sealed urine container

- ◇ Make sure that the sterile seal is intact before collecting. (Picture 1)
- ◇ Break the seal right before giving the container to the patient for urine collection.
- ◇ Within 4 minutes after the void, check the temperature is within the acceptable range (32°C to 38°C). Any urine specimen found outside the range, please recollect again.



- ◇ Please collect **2 urine specimens** with at least 20 mL each.
Note: In case of presumptive positive, 1 specimen will be used to send for confirmatory testing upon request by in-charge clinician.
- ◇ After urine collection and temperature check, the collector will then seal the container in the manner as shown in Picture 2. Then both the collector and patient must date, time and initial the seal.
- ◇ Send the sealed container to the laboratory.
- ◇ The laboratory staff will inspect the seal if it is kept intact during transportation. Any signs of tampering will be rejected by the laboratory.

Refer to changes highlighted yellow

STOOL COLLECTION

- ◇ Urinate before collecting the stool to avoid any urine in the stool sample. Do not urinate while passing the stool.
- ◇ Put on gloves before handling stool.
- ◇ Pass stool (but no urine) into a dry container. A plastic basin may be given that is able to be placed under the toilet seat to catch the stool.
- ◇ Either solid or liquid stool can be collected.
- ◇ Do not collect the sample from the toilet bowl.
- ◇ Do not mix toilet paper, water, or soap with the sample.
- ◇ Fill each container at least 1/3 full. **Do not overfill.**
- ◇ Wash hands thoroughly after removing glove.

SPUTUM COLLECTION

The specimen should be evaluated and recollected if sample appears to be saliva.


- ◇ Early morning specimens are preferred.
- ◇ Obtain the specimen when the cough is productive.
- ◇ Obtain sputum specimens before antibiotic therapy is initiated.
- ◇ Do not obtain a sputum specimen immediately after a patient has eaten, chewed gum, smoked, brushed teeth or uses mouthwash.
- ◇ Have patient rinse mouth with water to remove gross saliva before attempting cough. (Do not rinse with tap water if Legionella or AFB cultures are ordered)
- ◇ Send all specimens to the laboratory immediately in sterile screw-cap containers.

Mycobacterium tuberculosis (TB) Sputum Culture:

- ◇ Submit three first-morning expectorate sputum samples, collected on three different days. Label the sample accordingly, i.e. "1st Sample (Date: dd/mm/yy)", "2nd Sample (Date: dd/mm/yy)" and "3rd Sample (Date: dd/mm/yy)".
- ◇ At least 3 mL sputum is required for AFB cultures.
- ◇ Refrigerate the samples


Refer to changes highlighted yellow

NASOPHARYNGEAL SWAB COLLECTION



- ◇ Slightly tilt patient's head.
- ◇ Insert flexible swab through nares parallel to palate until:
 - a. Resistance is met, OR
 - b. Distance is equivalent to half the distance from the patient's ear to their nostril
- ◇ Gently rub and roll the swab.
- ◇ Leave swab in place for several seconds to absorb secretions.
- ◇ Slowly remove the swab while rotating it and immediately place in sterile tube (see individual test for specimen container requirement).

NASAL SWAB COLLECTION



- ◇ Slightly tilt patient's head.
- ◇ Insert swab about 2cm into a nostril.
- ◇ Gently rotate swab against the nasal wall for 5-10 times. Repeat collection procedure with second nostril.
- ◇ Slowly remove the swab and immediately place in sterile tube (see individual test for specimen container requirement).

SPECIMEN COLLECTION FOR FUNGUS

- ◇ **Skin:**
Wipe area with 70% alcohol and allow to dry completely. Scrap skin area using a scalpel blade gently to collect infected scales into a sterile container.
- ◇ **Scalp and Hair:**
Wipe area with 70% alcohol and allow to dry completely. Scrap scalp area gently to collect infected material (if scalp) or pull hair from affected scalp lesion, and place into a sterile container.
- ◇ **Nails:**
Wipe top and underside of nail with 70% alcohol and allow to dry completely. Collect nail clippings and place into a sterile container.

Other specimens may be submitted in a sterile container, tube or envelope; do not add saline or any other fluid to the container.

Submit specimens within 24 hours. Store and transport specimens at room temperature only.

SURGICAL PATHOLOGY TISSUE COLLECTION

- ◇ Specimens should be immersed in 10% formalin as soon as possible but within 1 hour of the biopsy or resection procedure. Use at least 3-4 times formalin to tissue volume.
- ◇ The following information must be legibly recorded on the label sticker on the specimen container (not the lid):
 - Patient's full name
 - Medical Record Number or other unique identifier
 - Date when specimen was obtained
 - Name of physician
 - Specimen type (and anatomical origin of site)

NOTE: HISTOPATHOLOGY SPECIMEN ARE OUTSOURCED TO A REFERRAL LAB

SUREPATH PAP TEST COLLECTION**Broom-Like Device Protocol**

1. Obtain an adequate sampling from the cervix using the Rovers Cervex-Brush. Insert the central bristles of the broom into the endocervical canal deep enough to allow the shorter bristles to fully contact the ectocervix. Push gently, and rotate the broom in a clockwise direction five times.
2. Put the broom as quickly as possible into the BD Sure Path[®] solution vial by dropping the detachable head of the device into the vial. Discard the collection device.
3. Tighten the cap so that the torque line on the cap passes the torque line on the vial.
4. Record the patient's name and ID number on the vial, and the patient information and medical history on the cytology requisition form. Send the SurePath vial to the lab for processing.

NOTE: SUREPATH SPECIMEN ARE OUTSOURCED TO A REFERRAL LAB

Refer to changes highlighted yellow


SPECIMEN LABELLING

◇ The person collecting the blood sample is responsible for transcribing the following info or at least two identifiers (patient’s name, and date of birth or PRN or Doc No.) on all sample containers:

- Patient’s name
- Patient’s MRN/PRN
- Patient’s document number (IC No. or passport)
- Date and year of birth
- Sex of patient
- Computer assigned lab # (if ordered through Vesalius)

Indicate source or site of sample if several samples are to be collected, e.g. swabs for culture, tissues

- ◇ Label tubes in the presence of the patient and right **after** the collection of the samples. Do not practice pre-labelling of the sample containers or tubes.
- ◇ Failure to properly label the tubes will require the specimens to be redrawn or recollection.
- ◇ The instructions below apply to all specimen tubes and containers (blood and non-blood):




Label the tube by placing the patient sticker straight vertically

Leave a window when labelling, to view the volume and quality of sample.


Cut sticker if the size covers the tube (must leave at least HIS-generated barcode, MRN, Full Name and Doc No.)

Stick patient sticker in the **same orientation** over the existing container label (e.g. from the top to bottom for blood tubes)


Do NOT label tube improperly as shown below




Twisted




Horizontal




Crinkled



Bottom-up (blood tube)



Lower part of barcode covered during testing, analyser unable to read



SPECIMEN HANDLING AND TRANSPORT

- ◇ All **IN HOUSE TEST** description and requirements related to type of blood collection devices, preservatives, sample volumes, etc. of different laboratory tests can be found in the test catalogue.
- ◇ For any special requests / send out test, or if in doubt for what container to use, please contact or phone the laboratory (office hours) to confirm type of blood container, blood volume, preservative and transport requirements. Microtainer systems are not ideal for send out test, it is available for pediatric and neonatal patient only.
- ◇ All specimens must be in properly sterile closed/capped container to avoid leakage. All specimen must be handled with care and treated as potentially infectious.

TRANSPORTATION OF SPECIMEN

- ◇ Send specimens to the laboratory promptly to maintain their integrity. Unspun or unprocessed (uncentrifuged) blood specimens must be delivered **within 2 hours of collection**. Failure to separate red cells from serum or plasma within 2 hours may result in inaccurate test results.
- ◇ Do **not** refrigerate unspun blood specimens before transport, as this can lead to elevated potassium levels.
- ◇ For specimens requiring protection from light, tightly wrap aluminum foil or a paper towel around the tube (but not the stopper). Ensure the tube is shielded from light during storage and transport.
- ◇ All specimens **MUST** be placed in sealed biohazard bag(s). Separate different sample types in different biohazard bag.
- ◇ Place specimen in the Ziploc portion of the specimen bag. Completed request form is to be placed in the outside pocket.
- ◇ Special requirement for transport can be found in the test catalogue for individual test (e.g. specimen requires chilling for Ammonia test).
- ◇ If specimen does not have a specific storage requirements, please store at **room temperature** before courier pick-up.
- ◇ All specimen transported to the laboratory by hand or courier must be transported in sealed biohazard, leak-proof, puncture resistant container tightly closed before transportation.
- ◇ All specimen transported to the laboratory by pneumatic tube must be cushioned with bubble wrap or any suitable padding inside the carrier capsule. **Do NOT send blood gas, surgical tissue, body fluids,**

Refer to changes highlighted yellow

bone marrow specimens and amniotic fluid for cytogenetic examination via Pneumatic Tube Systems.

REJECTION OF SPECIMEN

- ◇ JPMC laboratory will reject unacceptable specimen and inform the respective wards or clinics to request for a recollection for the following:
 - Mislabeled and Unlabeled specimens
 - Insufficient amount for the test requested
 - Haemolysed specimens
 - Inappropriate specimen container or transport media
 - Anticoagulated specimens containing blood clot
 - Improper storage or transport requirement
 - Discrepancies between requisition form and labeled tube
 - Specimen collected at the wrong time (e.g. drug level)
 - Specimen contamination or evidence of leakage

PRECIOUS SPECIMEN

- ◇ Special consideration will be made for precious specimens, where re-collection is not possible, such as normally Sterile Body Fluids (i.e.: pericardial, peritoneal, CSF), Bone Marrow, Fine Needle Aspirates (FNA), Tissue Biopsies/Surgical Specimens, and special timed specimen

The laboratory will proceed with testing only when the ordering physician agrees to accept responsibility for the specimen's identity, signs the Laboratory Precious Specimen Form and returns the form to the laboratory. This information will be recorded and included in the laboratory report.

LABORATORY REPORTING

- ◇ Laboratory results are treated with strict confidentiality. Access to laboratory results is restricted by password use and can only be viewed by the ordering location.
- ◇ Patient test results are available on-line after validation by the laboratory. In case of IT downtime, patient test results will be reported on paper and dispatched to the ordering location.

AMENDED REPORT

- ◇ Report that has been revised after release of results will be immediately informed to the requesting doctor or on-call doctor or primary nurse-in-charge.
- ◇ Revised result will be indicated on the report with a remark.
- ◇ The clinic or ward is responsible in ensuring that the original report is discarded, and the latest revised report is attached into the patient's medical record.

OUTSOURCED LABORATORY TESTS

- ◇ For tests that are not available in-house, specimen will be outsourced to an accredited referral laboratory.
- ◇ Laboratory report from referral laboratory will be available on-line through the Laboratory Information System and original hardcopy report will be dispatched to the ordering location.
- ◇ For outsourced tests turn around time, please contact Central Specimen Receiving and Management section at extension 2127.

CRITICAL RESULTS

- ◇ Critical result will be communicated immediately by phone to the requesting doctor. If requesting doctor is not available, on-call doctor or primary nurse-in-charge will be notified according to hospital policy.

Refer to changes highlighted yellow

CRITICAL RESULT LIST

BIOCHEMISTRY & IMMUNOLOGY CRITICAL RESULT			
Test	Low	High	Unit
Amikacin	-	> 8 (Trough) > 35 (Peak)	mg/L
Bilirubin, total (>1 month to < 1 year)	-	≥ 257	umol/L
Calcium	< 1.75	> 3.25	mmol/L
Creatine kinase	-	≥ 10000	U/L
Cyclosporine	< 100	> 800	ng/mL
Digoxin	-	> 2	ng/mL
Gentamicin	-	> 2 (Trough) > 12 (Peak)	mg/L
Glucose	< 2.5	> 25.0	mmol/L
Glucose (0 day to 2 months)	< 2.6	> 7.3	mmol/L
Glucose (2 months to 12 years)	< 3.3	> 25.0	mmol/L
Lactate	-	> 5	mmol/L
Magnesium	< 0.4	-	mmol/L
Osmolality (Serum)	≤ 190	≥ 390	mOsm/kg
Phenytoin	-	> 20	mg/L
Potassium	< 2.5	> 6.0	mmol/L
Sodium	< 120	> 160	mmol/L
Tacrolimus	< 4	> 12	ng/mL
Troponin-I	-	> 15.6 (Female) > 34.2 (Male)	ng/L
Valproic acid	-	> 100	mg/L
Vancomycin	-	> 10 (Trough) > 80 (Peak)	mg/L
HAEMATOLOGY CRITICAL RESULT			
Test	Low	High	Unit
APTT	-	> 100	sec
Fibrinogen	< 80	-	mg/dL
PT-INR	-	> 5.0	-
WBC Count	< 2.0	> 50.0	x 10 ⁹ /L
Absolute Neutrophil Count	< 0.5	-	x 10 ⁹ /L
HB (0-7 weeks)	≤ 60	> 240	g/L
HB (> 7 weeks)	≤ 60	≥ 200	g/L
Platelet Count	< 50	> 800	x 10 ⁹ /L
Malaria Parasites	Positive		
Blast	Present (1st presentation or at relapse)		
MICROBIOLOGY CRITICAL RESULT			
Positive for blood culture, AFB smear and CSF gram stain			

BLOOD TRANSFUSION POLICY

CROSSMATCH REQUEST

Specimen: Plasma (EDTA, Minimum volume: 1 mL)
 Collection: Freshly drawn sample is required
Person collecting the blood must sign/initial and date on the patient sticker on both blood sample and request form
 Specimen Stability: 3 days at 2-8°C. If plasma is separated, then 7 days at 2-8°C.
 Rejection Criteria: Grossly haemolysed
 Method: Column Agglutination using Ortho Vision System or Semi-Automated Method
 Useful for: Performed prior to blood transfusion to determine compatibility between donor’s blood and blood of the recipient
 Performed: Daily
 TAT: 1 hour, STAT – 40 mins, Emergency Release – 10 mins
 Section: Blood Transfusion (ext 2320)

Type in Blood Product Code, followed by number of products required.	
Blood product	Blood product Code
Whole blood	WB
Packed cells	PC
Leukodepleted Packed cells	LDPC
Fresh Frozen Plasma	FFP
Cryoprecipitate	CRYO
Platelets	PLTC
Plateletpheresis	APH

- ◇ To order crossmatch in the electronic system, please insert blood product code and no. of units required. For example:
 - To request for 4 units of packed cells, the code is PC4.
 - To request for 1 unit of Plateletpheresis, the code is APH1.

- ◇ All blood product requisitions must be ordered and signed by physicians using the BLOOD BANK REQUEST FORM.

- ◇ If ward/clinic is using the Vesalius system to request test, they will still need to fill in the BLOOD BANK REQUEST FORM in addition to the electronic request and ensure that the Vesalius patient sticker with the lab number is pasted on the blood bank request form.

Refer to changes highlighted yellow

INVESTIGATION OF ADVERSE TRANSFUSION REACTION

Specimen:	<u>Post-transfusion:</u> <ul style="list-style-type: none"> - Plasma (EDTA, Minimum volume: 1 mL) - Next voided Urine - Donor's blood bag returned in a biohazard bag - Transfusion Reaction Notification Form - Plain/SST Tube (IF requested by the physician for Bilirubin test) - Blood Culture bottles (Only IF doctor suspects bacterial contamination of blood product).
Rejection Criteria:	Haemolysed specimen (a recollection of post-transfusion blood specimen is required to ensure haemolysis is not due to collection procedure)
Method:	Perform initial, primary and/or secondary investigation with pre- and post-transfusion patient's sample
Performed:	Daily
Analytical Time:	STAT

◇ The attending physician or nurse in-charge must immediately inform the Blood Transfusion Section if any transfusion reaction occurs, and transfusion of the blood product must be discontinued if major symptoms are observed. Reconfirm the identification of patient and the blood product infused with the Blood Bank.

◇ Transfusion reactions include:

- Fever without chills (associated with Haemolytic Transfusion Reaction)
- Skin symptoms like hives (urticaria) or itching
- Pain (chest, abdomen and/or flank)
- Acute hypotension or acute hypertension
- Respiratory diseases (dyspnea, tachypnea, wheezing or hypoxemia)
- Nausea and/or vomiting
- Darkened urine or jaundice
- Flushing
- Bleeding
- Anaphylaxis

◇ Obtain post-transfusion samples from the patient immediately (or as soon as possible in the case of urine samples) for workup. Complete a Transfusion Reaction Notification Form with required information. A copy of the form must accompany the samples to the Blood Bank.

Note: Be careful to minimize mechanical haemolysis when taking blood samples. Tubes must also be

Refer to changes highlighted yellow

properly labelled.

- ◇ Send blood samples, Crossmatch / Component Transfusion form, Transfusion Reaction Notification Form and blood bag with infusion set and tag attached to the Blood Bank.
- ◇ Do not transfuse any more products until the reaction investigation is completed.

Delayed Transfusion Reaction

- ◇ Haemolytic transfusion reactions may not be apparent for several days following transfusion. Most common signs are fever, unexpected fall or less than expected rise in haemoglobin, and jaundice.
- ◇ If any delayed haemolytic transfusion reactions are observed by the nurses, notify the attending physician immediately for any medical treatment necessary for the patient.
- ◇ A Transfusion Reaction Notification Form must be submitted again to the Blood Bank together with the patient's post-transfusion specimens (1 x EDTA, 50 mL urine).

EMERGENCY RELEASE OF BLOOD UNITS

Specimen:	Plasma (EDTA, Minimum volume: 1 mL)
Rejection Criteria:	Grossly haemolysed
Method:	Tube method and Ortho Vision system AFTER blood products have been issued
Useful for:	For the immediate release of un-crossmatched group and type specific blood products or group O blood products, either Rh-positive or Rh-negative

- ◇ The attending physician must issue a written order using Emergency Release of Blood Form specifying blood group (if known) and number of units required. The form must be completely prescribed and signed by the attending physician.
- ◇ At the same time, the attending staff nurse **must** inform the laboratory of the urgency of blood product/s needed to avoid unnecessary delays.
- ◇ One EDTA tube of blood must be drawn from the patient and sent down to the laboratory together with all the required forms.
- ◇ Laboratory staff will inform the nurse in-charge of available units and issue blood products as requested and will continue with the crossmatch procedure simultaneously.

Refer to changes highlighted yellow

- ◇ Should there be any incompatibilities noted with the crossmatch procedure, laboratory staff will immediately notify the nurse in-charge, and the transfusion process must be stopped. The blood unit must then be returned to the laboratory as soon as possible for adverse transfusion reaction investigation.

Note:

- There is no emergency release for FFP or Cryoprecipitate products as thawing will take up to 20 minutes from the receipt of the request.
- Platelet concentrate stock depends on availability of platelet donors. Pure platelets are harvested via Apheresis process, which may take more than a day to obtain.

ANNEX A. TEST CATALOGUE (IN-HOUSE)

25-OH VITAMIN D TOTAL, Code: OHD

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) or Plasma (EDTA, purple top, 4mL)

Collection: NA

Specimen Stability: 12 days at 2-8°C

Rejection Criteria: NA

Method: Chemiluminescent Microparticle Immunoassay (CMIA)

Useful for: Diagnosis of Vitamin D deficiency and Hypervitaminosis D

Reference range: See laboratory report

Performed: Every Tuesday and Friday

TAT: 3 days

Section: Immunology (ext 2407)

ABSCESS SWAB CULTURE & SENSITIVITY, Code: ABS

Specimen: Abscess (Gel Swab/Sterile container/Syringe with needle removed)

Collection: Indicate source of specimen

Specimen Stability: Room temperature, 24 hours

Rejection Criteria: Dry swab, samples with formalin

Method: Conventional culture

Useful for: Diagnosis of bacterial infection

Reference range: -

Performed: Daily

TAT: 48-72 hours

Section: Microbiology (ext 2436)

ACINETOBACTER SCREENING, Code: Refer to note below

OTHER TEST NAME: STERILITY TESTING

Specimen: Any (Gel swab) / Sterile container (for urine and sputum)

Collection: Indicate source of specimen

Specimen Stability: Room temperature, 24 hours

Rejection Criteria: Dry swab, Source of specimen not indicated if more than one swab

Method: Conventional culture

Useful for: Detection of *Acinetobacter baumannii*

Reference range: -

Performed: Daily

TAT: 48-72 hours

Section: Microbiology (ext 2436)

NOTE: Code for sterility test is STER (no. of specimen)
(i.e. For one specimen, code is STER1. For two specimens, code is STER2. And so on and so forth).

ACTIVATED PARTIAL THROMBOPLASTIN TIME, Code: APTT

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line indicated on the tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Coagulometric Measurement
 Useful for: Monitoring heparin therapy and screening test for clotting factors
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

ACTIVATED PARTIAL THROMBOPLASTIN TIME RATIO, Code: APTRATIO

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line indicated on the tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Coagulometric Measurement
 Useful for: NA
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

APTT 50% CORRECTION, Code: APT50

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line indicated on the tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Coagulometric Measurement
 Useful for: Follow up test to investigation cause of unexplained prolong APTT. Differentiates between possibility of factor deficiency or acquired coagulation inhibitors.
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day
 Section: Haematology (ext 2408)

Refer to changes highlighted yellow

ALANINE AMINOTRANSFERASE (SGPT), Code: ALT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: IFCC, NADH without pyridoxal phosphate
 Useful for: Aid in diagnosis and treatment of certain liver diseases (e.g viral hepatitis and cirrhosis)
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

ALBUMIN, Code: ALB

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) or Plasma (EDTA, purple top, 4mL)
 Collection: NA
 Specimen Stability: 7 days at room temperature and at 2-8°C
 Rejection Criteria: NA
 Method: Colorimetric with Bromcresol Green
 Useful for: Aid in the diagnosis and treatment of numerous diseases involving primarily in liver or kidneys.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

ALBUMIN & CREATININE RATIO, Code: ACR

Specimen: Urine (Sterile container or 24hr urine container)
 Collection: No preservative needed for random urine.
 Preserve urine using 20 to 30 mL of 6M HCL for 24 hours urine and store refrigerated during collection period
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Renal function test
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

ALKALINE PHOSPHATASE, Code: ALP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: IFCC, Para-nitrophenyl Phosphate
 Useful for: **Aid in the diagnosis and treatment of liver, bone, parathyroid and intestinal diseases.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

ALPHA-FETOPROTEIN, Code: AFP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: A tumour marker for hepatocellular carcinoma and testicular cancer.
Aid in detection of fetal open neural tube defects (NTD).
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 2 days
 Section: Immunology (ext 2407)

AMIKACIN, Code: AMIK

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

AMIKACIN LEVEL (PEAK) , Code: AMIKP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Collect 30 min after end of IV infusion or 1 hour after IM injection
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

AMIKACIN LEVEL (TROUGH) , Code: AMIKT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Collect specimen immediately before next dose
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

AMMONIA, Code: AMON

Specimen: Plasma (**EDTA** purple top, 4 mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: **Place on ice, send to the lab immediately within 15 minutes for rapid centrifugation**
***(please inform lab at least 1 hour prior to collection)**
 Specimen Stability: Assay immediately
 Rejection Criteria: **Specimen is not placed on ice**
 Method: Enzymatic with Glutamate Dehydrogenase
 Useful for: **Liver function assessment.**
Screening test for amino acid disorders, organic acid disorders and urea cycle disorders.
 Reference range: See laboratory report
 Performed: Daily
 TAT: **2 hours**
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

AMYLASE, Code: AMY

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 24 hours at room temperature & 2-8°C
 Rejection Criteria: **NA**
 Method: IFCC, Enzymatic with CNPG3 Substrate
 Useful for: **Diagnosis and treatment of pancreatitis**
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Biochemistry (ext 2407)

AMYLASE (URINE RANDOM), Code: UAM

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 24 hours at room temperature, 3 days at 2-8°C
 Rejection Criteria: Urine collected with acidic preservatives.
 Method: Enzymatic/Colorimetric with EPS
 Useful for: **Diagnosis and treatment of pancreatitis**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

AMYLASE (24 HOURS URINE), Code: UAM24

Specimen: Urine (24hr urine container)
 Collection: Timed or 24 hours specimens with no preservatives
 Specimen Stability: 24 hours at room temperature, 3 days at 2-8°C
 Rejection Criteria: Urine collected with acidic preservatives.
 Method: Enzymatic/Colorimetric with EPS
 Useful for: **Diagnosis and treatment of pancreatitis**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

ANTIBODY SCREENING, Code: BAB

OTHER TEST NAME: INDIRECT ANTIGLOBULIN TEST

Specimen: Plasma from Whole Blood (EDTA, **Minimum volume: 1 mL**)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: Grossly Haemolysed
 Method: Column Agglutination using Ortho Vision System
 Useful for: To detect as many clinically significant antibodies as possible
 Performed: Daily
 TAT: 1 day
 Section: Blood Transfusion (ext 2320)

ANTI-CCP (Cyclic Citrulinated Peptide) , Code: CCP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) **or Plasma (EDTA, purple top, 4mL)**
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Aids in diagnosis of Rheumatoid Arthritis
 Reference range: See laboratory report
 Performed: Monday, Wednesday & Friday
 TAT: 2 days
 Section: Immunology (ext 2407)

ANTI-THYROGLOBULIN , Code: ATHY

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) **or Plasma (EDTA, purple top, 4mL)**
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Aids in diagnosis of autoimmune thyroid diseases, Hashimoto's Thyroiditis, Primary Myxedema and Graves' disease. Useful in follow-up of patients with thyroid carcinomas.
 Reference range: See laboratory report
 Performed: Wednesday & Saturday
 TAT: 2-3 days
 Section: Immunology (ext 2407)

ANTI-THYROID PEROXIDASE , Code: TPO	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) or Plasma (EDTA, purple top, 4mL)
Collection:	NA
Specimen Stability:	3 days at 2-8°C
Rejection Criteria:	NA
Method:	Chemiluminescent Microparticle Immunoassay (CMIA)
Useful for:	Aids in diagnosis of autoimmune thyroid diseases, Hashimoto's Thyroiditis, Primary Myxedema and Graves' disease.
Reference range:	See laboratory report
Performed:	Wednesday & Saturday
TAT:	2-3 days
Section:	Immunology (ext 2407)
ANTI-THYROID RECEPTOR ANTIBODY , Code: TRAB	
Specimen:	Serum (SST, yellow top, 5mL)
Collection:	NA
Specimen Stability:	3 days at 2-8°C
Rejection Criteria:	NA
Method:	Chemiluminescent Microparticle Immunoassay (CMIA)
Useful for:	Aids in differential diagnosis of Graves' disease. Useful in monitoring anti-thyroid therapy and predicting remission.
Reference range:	See laboratory report
Performed:	Wednesday & Saturday
TAT:	2-3 days
Section:	Immunology (ext 2407)
APT-DOWNEY TEST, Code: APT	
Specimen:	Stool (Sterile container with attached spatula) or vomit (Sterile container)
Collection:	Fresh with minimal debris
Specimen Stability:	Test immediately
Rejection Criteria:	NA
Method:	APT-Downey test
Useful for:	To distinguish between maternal (adult type) and infant's (fetal type) hemoglobin in a grossly bloody stool.
Reference range:	N/A
Performed:	Daily
TAT:	1 day
Section:	Biochemistry (ext 2407)

Refer to changes highlighted yellow

ARTERIAL/VENOUS BLOOD GAS, Code: GAS (for arterial), VBG (for venous)	
Specimen:	Whole blood (Heparinized syringe, 1 mL)
Collection:	Do not expose specimen to air and mix well to prevent clotting. Chill specimen on ice and send immediately to lab. Do not send through pneumatic tube.
Specimen Stability:	Test immediately
Rejection Criteria:	The specimen was neither placed on ice nor sent immediately, evidence of clotting, presence of air bubbles and specimen sent through pneumatic tube.
Method:	Abbott i-STAT
Useful for:	Aid in diagnosis, monitoring and treatment of respiratory and metabolic acid-base disturbances.
Reference range:	See laboratory report
Performed:	Daily
TAT:	STAT
Section:	Biochemistry (ext 2407)
ASPARTATE AMINOTRANSFERASE (SGOT), Code: AST	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
Collection:	NA
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	IFCC, Enzymatic without pyridoxal phosphate
Useful for:	Liver profile assessment. Diagnosis of Acute Myocardial Infarct.
Reference range:	See laboratory report
Performed:	Daily
TAT:	1 day
Section:	Biochemistry (ext 2407)
BENCE JONES PROTEIN, Code: BJP	
Specimen:	Urine (Sterile container)
Collection:	Fresh morning urine
Specimen Stability:	Test immediately
Rejection Criteria:	NA
Method:	Heat test
Useful for:	Screening test for multiple myeloma and amyloidosis
Reference range:	Not detected in normal individuals. Plasma cells leukemia will need to do immunoelectrophoresis technique because it is negative for Bence Jones protein
Performed:	Office hours only
TAT:	1 day
Section:	Biochemistry (ext 2407)

Refer to changes highlighted yellow

BICARBONATE, Code: HCO3

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Not recommended as add-on as a consequent decrease in the CO2 value of up to 6 mmol/L can occur in an hour once the specimen has been exposed to ambient air.
 Specimen Stability: 2 hours at room temperature (Tightly capped)
 Rejection Criteria: NA
 Method: Enzymatic with PEP Carboxylase
 Useful for: Evaluate acid-base imbalances
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

BILIRUBIN-DIRECT, Code: BILD

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Diazo reaction
 Useful for: Differential diagnosis of jaundice
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

BILIRUBIN-TOTAL, Code: BILT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Diazonium salt
 Useful for: Assessment of liver function and measurement of unconjugated bilirubin.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

BLOOD CULTURE (Aerobic or Anaerobic), Code: BC

Specimen: Whole blood (Bactec blood culture bottles)
 Collection: 8 -10 mL of blood. Recommended to draw two blood culture sets per septic episode from different sites at an interval of 45 mins – 1 hour between each set
 Specimen Stability: Room temperature for 4 hours
 Rejection Criteria: Keep at room temperature for more than 4 hours
 Method: Automated Bactec (fluorescence) & Conventional culture
 Useful for: Diagnosis of septicaemia
 Reference range: No growth after 5 days incubation
 Performed: Daily
 TAT: 5 days
 Section: Microbiology (ext 2436)

BLOOD CULTURE for PAEDIATRIC (Aerobic) Code: BCP

Specimen: Whole blood (1-3 mL, Bactec blood culture bottle)
 Collection: NA
 Specimen Stability: Room temperature for 4 hours
 Rejection Criteria: Keep at room temperature for more than 4 hours
 Method: Automated Bactec (fluorescence) & Conventional culture
 Useful for: Diagnosis of septicaemia
 Reference range: No growth after 5 days incubation
 Performed: Daily
 TAT: 5 days
 Section: Microbiology (ext 2436)

BLOOD CULTURE for YEAST IDENTIFICATION and SENSITIVITY, Code: BCY

Specimen: Whole blood (Bactec Myco/F Lytic blood culture bottle)
 Collection: 1 – 5 mL of blood
 Specimen Stability: Room temperature for 4 hours
 Rejection Criteria: Keep at room temperature for more than 4 hours
 Method: Automated Bactec (fluorescence) & Conventional culture
 Useful for: Diagnosis of yeast infection in the blood
 Reference range: No growth after 30 days incubation
 Performed: Daily
 TAT: 30 days
 Section: Microbiology (ext 2436)

BLOOD FILM, Code: BF

Specimen: Whole blood (EDTA purple top)
 Collection: 4ml for EDTA tube and 500ul for microtainer
 Specimen Stability: 12 hours from the time of collection
 Rejection Criteria: Haemolysed and/or clotted
 Method: Microscopy
 Useful for: To provide information of the morphology of blood cells
 Reference range: See Laboratory Report
 Performed: Office hours only
 TAT: 2-3 days
 Section: Haematology (ext 2408)

BLOOD GROUPING, Code: BG

Specimen: Whole Blood (EDTA purple top, Minimum volume: 1 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: Grossly Haemolysed
 Method: Column Agglutination using Ortho Vision System or Manual Tube Method
 Useful for: Blood group determination
 Performed: Daily
 TAT: 1 day
 Section: Blood Transfusion (ext 2320)

CA 125, Code: CA1

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Tumour marker for ovarian cancer**
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 2 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

CA 15-3, Code: CA5

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Tumour marker for Stage II and III breast cancer**
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 2 days
 Section: Immunology (ext 2407)

CA 19-9, Code: CA9

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Tumour marker for pancreatic and hepatobiliary cancer**
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 2 days
 Section: Immunology (ext 2407)

CALCIUM TOTAL, Code: CA

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 3 weeks at 2-8°C
 Rejection Criteria: **NA**
 Method: Colorimetric method with Arsenazo III
 Useful for: Evaluation of calcium metabolism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CALCIUM IONISED, Code: CAI

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Send sample immediately to lab. Not recommended for add-on as exposure of sample to air will result in decrease of ionized calcium
 Specimen Stability: Test immediately
 Rejection Criteria: NA
 Method: Ion-selective electrode potentiometric
 Useful for: Evaluation of calcium metabolism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CALCIUM (URINE RANDOM), Code: UCA

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 4 days at 2-8°C
 Rejection Criteria: NA
 Method: Colorimetric method with Arsenazo III
 Useful for: Evaluation of calcium metabolism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CALCIUM & CREATININE RATIO (URINE RANDOM), Code: RCCR

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 4 days at 2-8°C
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Screening tool for hypercalciuria
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CALCIUM (24 HOURS URINE), Code: UCA24

Specimen: Urine (24hr urine container)
 Collection: Preserve urine using 20 to 30 mL of 6M HCL and store refrigerated during collection period
 Specimen Stability: 4 days at 2-8°C
 Rejection Criteria: NA
 Method: Colorimetric method with Arsenazo III
 Useful for: Evaluation of calcium metabolism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CALCIUM & CREATININE RATIO (24 HOURS URINE), Code: CCR

Specimen: Urine (24hr urine container)
 Collection: Preserve urine using 20 to 30 mL of 6M HCL and store refrigerated during collection period
 Specimen Stability: 4 days at 2-8°C
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Screening tool for hypercalciuria
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CARBAMAZEPINE, Code: CARB

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CARCINOEMBRYONIC ANTIGEN, Code: CEA

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Tumour marker in colon and pancreatic cancer**
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 2 days
 Section: Immunology (ext 2407)

CARDIAC ENZYMES, Code: CAR

Tests include Creatine Kinase, Lactate Dehydrogenase and Aspartate Aminotransferase

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: **3 days at 2-8°C**
 Rejection Criteria: **NA**
 Method: Refer to individual test method
 Useful for: **Assessment of cardiac disorders**
 Reference range: See laboratory report
 Performed: **Daily**
 TAT: **1 day**
 Section: Biochemistry (ext 2407)

CHLAMYDIA TRACHOMATIS/NEISSERIA GONORRHOEAE DNA, CODE: CTNG

Specimen: Endocervical Swab or First Void Urine (7 mL)
 (Special swab collection kit from the Laboratory)
 Collection: NA
 Specimen Stability: 4 days at 2-8°C (Urine)
 Rejection Criteria: Swab not in Xpert transport media
 Method: NAAT
 Useful for: To aid in the diagnosis of chlamydial and gonorrhoeal disease
 Reference range: Not Detected
 Performed: Office hours only
 TAT: 1-3 days
 Section: Virology and Serology (ext 2322)

Refer to changes highlighted yellow

CHLORIDE, Code: CL

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: **Assessment of electrolyte imbalance**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

CHLORIDE (URINE RANDOM), Code: UCL

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation / assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CHLORIDE (24 HOURS URINE), Code: UCL24

Specimen: Urine (24hr urine container)
 Collection: No preservative needed. Store refrigerated during collection period
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation / assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CHOLESTEROL TOTAL, Code: CHO

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: N/A
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Enzymatic
 Useful for: Diagnosis of hyperlipidemia and cardiovascular risk assessment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CLOSTRIDIUM DIFFICILE TOXIN, Code: CDST

Specimen: Stool (Sterile container)
 Collection: -
 Specimen Stability: Within 24 hours (Room Temperature), 72 hours at 2-8°C
 Rejection Criteria: Specimen exceeds stability, specimen kept in formalin/fixative
 Method: Rapid immunochromatographic assay
 Useful for: Detects C.difficile antigen and toxin
 Reference range: Negative
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

COMPLEMENT C3 , Code: C3

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Immunospectrophotometric
 Useful for: Aids in determination of inherited or acquired deficiencies, and diagnosis of inflammatory and necrotic disorders.
 Reference range: See laboratory report
 Performed: Monday, Wednesday & Friday
 TAT: 2 days
 Section: Biochemistry (ext 2407)

COMPLEMENT C4 , Code: C4

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4mL)
 Collection: NA
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: NA
 Method: Immunoturbidimetric
 Useful for: Aids in determination of inherited or acquired deficiencies, and diagnosis of inflammatory and necrotic disorders.
 Reference range: See laboratory report
 Performed: Monday, Wednesday & Friday
 TAT: 2 days
 Section: Biochemistry (ext 2407)

CORTISOL, Code: RCO

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Aids in diagnosis and treatment of adrenal disorders e.g Cushing's syndrome and Addison's disease
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days, STAT
 Section: Immunology (ext 2407)

CORTISOL (URINE RANDOM) , Code: UCO

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Screening test for Cushing's syndrome and Addison's disease
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days, STAT
 Section: Immunology (ext 2407)

CORTISOL (24 HOURS URINE), Code: UCO24

Specimen: Urine (24hr urine container)
 Collection: No preservative needed. Store refrigerated during collection period
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Screening test for Cushing’s syndrome and Addison’s disease
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days, STAT
 Section: Immunology (ext 2407)

C-REACTIVE PROTEIN, Code: CRP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 2 months at 2-8°C
 Rejection Criteria: NA
 Method: Turbidimetric/Immunturbidimetric
 Useful for: A diagnostic indicator of infections and inflammation and monitoring response to therapy.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CREATINE KINASE, Code: CK

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: IFCC with N-acetyl-L-cysteine
 Useful for: Assessment of skeletal & cardiac muscle disorders
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CREATINE KINASE-MB, Code: CMB

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Assist in diagnosis of myocardial infarction
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Immunology (ext 2407)

CREATININE, Code: CRT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Kinetic Alkaline Picrate
 Useful for: Aid in diagnosis and treatment of renal diseases and in monitoring renal dialysis
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

CREATININE (URINE RANDOM), Code: UCRT

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: NA
 Method: Kinetic Alkaline Picrate
 Useful for: Used to calculate creatinine clearance and as a calculation basis for measuring other urine analytes.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CREATININE (24 HOURS URINE), Code: UCRT24

Specimen: Urine (24hr urine container)
 Collection: **No preservative needed**
 Specimen Stability: **3 days at 2-8°C**
 Rejection Criteria: NA
 Method: Kinetic Alkaline Picrate
 Useful for: **Used to calculate creatinine clearance and as a calculation basis for measuring other urine analytes.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CREATININE CLEARANCE (24 HOURS URINE), Code: CC

Specimen: 1) Urine (24hr urine container)
 2) Serum (SST yellow top, 5 mL) or plasma (Li-Heparin green top, 4 mL)
 Note: Patient's weight (in kg) and height (in cm) are required for creatinine clearance calculation
 Collection: **No preservative needed**
 Specimen Stability: **3 days at 2-8°C**
 Rejection Criteria: **Only one specimen type received**
 Method: Calculated
 Useful for: Creatinine clearance evaluation
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

CRYPTOCOCCAL STAIN, Code: CPT

Specimen: CSF (Sterile container)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: -
 Method: Conventional
 Useful for: Diagnosis of meningitis
 Reference range: Negative
 Performed: Office hours only
 TAT: 1 day
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

CRYOGLOBULIN TEST, Code: CGB

Specimen: Serum (CAT, Red Top, 5mL)
 Collection: Must inform the laboratory a day before testing required.
 Incubated tubes and equipment must be used.
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: Icteric/Lipaemic
 Method: Incubation at 4°C for 72 hours and observed for precipitation
 Useful for: To help detect the presence of cryoglobulins in the blood to help determine or rule out potential causes of cryoglobulinemia
 Reference range: Negative
 Performed: Office hours only
 TAT: 3 days from time of collection
 Section: Blood Transfusion (ext 2320)

CSF ALBUMIN, Code: CSFALB

Specimen: Cerebrospinal Fluid (Sterile screw-capped container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: **NA**
 Method: Turbidimetric/Immunoturbidimetric
 Useful for: Assessment of CNS disease and infection
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Biochemistry (ext 2407)

CSF CULTURE only, Code: CSFC

Specimen: CSF (Sterile container)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: -
 Method: Light microscopy and conventional culture
 Useful for: Diagnosis of bacterial meningitis
 Reference range: Culture sterile
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

CSF FEME (GRAMSTAIN, MICROSCOPY, CELL COUNT), Code: CSFF

Specimen: CSF (Sterile container)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: Heavily blood-stained specimen
 Method: Light microscopy
 Useful for: Diagnosis of bacterial meningitis
 Reference range: Leukocytes 0-5 µl
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

CSF FEME & CULTURE, Code: CSFFC

Specimen: CSF (Sterile container)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: Heavily blood-stained specimen (for FEME)
 Method: Conventional culture
 Useful for: Diagnosis of bacterial meningitis
 Reference range: Culture sterile, leukocytes 0-5µl
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

CSF GLUCOSE, Code: CSFB

Note: Test is run together with CSF Protein

Specimen: Cerebrospinal Fluid (Sterile screw-capped container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: **NA**
 Method: Hexokinase/G-6-PDH
 Useful for: Assessment of CNS disease and infection
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

CSF PROTEIN, Code: CSFB

Note: Test is run together with CSF Glucose

Specimen: Cerebrospinal Fluid (Sterile screw-capped container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: **NA**
 Method: Turbidimetric with Benzethonium Chloride
 Useful for: Assessment of CNS disease and infection
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

CSF LACTATE, Code: CSFLAC

Specimen: Cerebrospinal Fluid (Sterile screw-capped container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: **NA**
 Method: Enzymatic, Lactic acid to pyruvate
 Useful for: Assessment of CNS disease and infection
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Biochemistry (ext 2407)

CSF LACTATE DEHYDROGENASE (LDH), Code: LDHCSF

Specimen: Cerebrospinal Fluid (Sterile screw-capped container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: **NA**
 Method: IFCC, Lactate to Pyruvate
 Useful for: Assessment of CNS disease and infection
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Biochemistry (ext 2407)

CYCLOSPORINE, Code: CYCLO

Specimen: Whole blood (EDTA, purple top, 3mL).
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Therapeutic Drug Monitoring for Organ Transplant
 Reference range: See laboratory report
 Performed: Daily, STAT
 TAT: 1 day
 Section: Immunology (ext 2407)

D-DIMER, Code: FDP

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line indicated on tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Immunological Measurement
 Useful for: Aid in the diagnosis of disseminated intravascular coagulation (DIC)
 Reference range: Refer Laboratory Report
 Performed: Daily, STAT
 TAT: 1 day
 Section: Haematology (ext 2408)

DEHYDROEPIANDROSTERONE SULPHATE, Code: DHEA

Specimen: Serum (SST yellow top, 5mL)
 Collection: NA
 Specimen Stability: 8 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Evaluation of androgen status
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

DENGUE SEROLOGY (NS1 Antigen + Antibody), Code: DEA

Specimen: Serum (SST yellow top, 5 ml)
 Collection: NA
 Specimen Stability: 2 weeks at 2-8°C
 Rejection Criteria: NA
 Method: Rapid immunoassay
 Useful for: Diagnosis of acute/past dengue infection
 Reference range: Negative
 Performed: Office hours only
 TAT: 1 day
 Section: Virology and Serology (ext 2322)

DIABETIC SCREEN, Code: DIA

Tests include fasting glucose and glycated haemoglobin (HbA1c)

Specimen: 1) Plasma (Sodium fluoride, grey top, 3mL) for Fasting Glucose
 2) Whole blood (EDTA, purple top, 3mL) for HbA1c
 Collection: Fasting for at least 8 hours
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: Clotted purple top. Not fasting or fasting for less than 8 hours.
 Method: Refer to individual test method
 Useful for: Diagnosis and monitoring of diabetes mellitus
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

DIABETIC STUDIES, Code: DST

Tests include electrolytes, urea, creatinine, fasting glucose, glycated haemoglobin (HbA1c) and urine microalbumin

Specimen: 1) Serum (SST, yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) for Electrolytes, Urea and Creatinine
 2) Plasma (Sodium fluoride, grey top, 3mL) for Fasting Glucose
 3) Whole blood (EDTA, purple top, 3mL) for HbA1c
 4) Urine (Sterile container) for Urine Microalbumin

Collection: Fasting for at least 8 hours

Specimen Stability: 5 days at 2-8°C

Rejection Criteria: Clotted purple top. Not fasting or fasting for less than 8 hours.

Method: Refer to individual test method

Useful for: Diagnosis and monitoring of diabetes mellitus

Reference range: See laboratory report

Performed: Daily

TAT: 1 day

Section: Biochemistry (ext 2407)

DIGOXIN, Code: DIG

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)

Collection: Draw blood 12 hours after oral administration or 6 hours after IV administration

Specimen Stability: 2 days at 2-8°C

Rejection Criteria: NA

Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)

Useful for: Therapeutic Drug Monitoring

Reference range: See laboratory report

Performed: Daily

TAT: 1 day, STAT

Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

DIRECT COOMB'S TEST, Code: DCT

OTHER TEST NAME: DIRECT ANTIGLOBULIN TEST

Specimen: Whole Blood (EDTA purple top, Minimum volume: 1 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: Grossly Haemolysed
 Method: Tube agglutination
 Useful for: Transfusion Reaction Investigation
 Diagnosis of haemolytic disease of the new born
 Diagnosis of autoimmune and drug induced haemolytic anemia
 Detection of passively acquired alloantibodies
 Reference Range: NA
 Performed: Daily
 TAT: 1 day
 Section: Blood Transfusion (ext 2320)

Du test, Code: DU

Specimen: Whole Blood (EDTA purple top, Minimum volume: 1 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: Grossly haemolysed
 Method: Tube method
 Useful for: Automatically performed to verify rhesus negative blood groups
 Reference Range: NA
 Performed: Daily
 TAT: 1 day
 Section: Blood Transfusion (ext 2320)

EAR SWAB CULTURE & SENSITIVITY, Code: RSCEAR

Specimen: Ear Swab (Gel swab)
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab
 Method: Conventional culture
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

ELECTROLYTES, Code: ELY

Tests include Sodium, Potassium, Chloride

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

ELECTROLYTES (URINE RANDOM), Code: ELY2

Tests include Urine sodium, Urine potassium, Urine chloride, Urine urea

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

ELECTROLYTES, CREATININE & UREA, Code: ECU

Tests include Sodium, Potassium, Chloride, Creatinine, Urea

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Assessment of renal function
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

ERYTHROCYTE SEDIMENTATION RATE, Code: ESR	
Specimen:	Whole blood (EDTA purple tube, 4 ml)
Collection:	NA
Specimen Stability:	4 hours at room temperature
Rejection Criteria:	Haemolysed and/or clotted, lipaemic
Method:	Photometric Rheology
Useful for:	An indicator of the presence and extent of inflammation and its response to treatment
Reference range:	Refer Laboratory Report
Performed:	Daily
TAT:	1 day
Section:	Haematology (ext 2408)
ESTRADIOL, Code: E2	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
Collection:	NA
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	Chemiluminescent Microparticle Immunoassay (CMIA)
Useful for:	Female: Evaluation of hypothalamic-pituitary-ovarian axis Male: Investigation of unexplained gynecomastia. Investigation of Infertility. Monitoring during IVF procedures to assess development of ovarian follicles
Reference range:	See laboratory report
Performed:	Every Monday and Thursday
TAT:	3 days, STAT
Section:	Immunology (ext 2407)
ETHANOL (BLOOD), Code: ETH	
OTHER TEST NAME: ALCOHOL (BLOOD)	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
Collection:	NA
Specimen Stability:	6 months at 2-8°C
Rejection Criteria:	NA
Method:	Alcohol Dehydrogenase
Useful for:	To detect presence and levels of alcohol
Reference range:	See laboratory report
Performed:	Daily
TAT:	1 day
Section:	Biochemistry (ext 2407)

Refer to changes highlighted yellow

ETHANOL (URINE), Code: UETH

OTHER TEST NAME: ALCOHOL (URINE)

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 30 days at 2-8°C
 Rejection Criteria: NA
 Method: Alcohol Dehydrogenase
 Useful for: To detect presence and levels of alcohol
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

EYE SWAB CULTURE & SENSITIVITY, Code: RSCEYE

Specimen: Eye Swab (Gel swab)
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab
 Method: Conventional culture
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

FAECAL OCCULT BLOOD, Code: OCI

Specimen: Stool (Sterile container with attached spatula)
 Collection: -
 Specimen Stability: 72 hours at 2-8°C
 Rejection Criteria: Specimens other than stool
 Method: Immunochromatographic Test
 Useful for: Detection of blood in stool
 Reference range: Negative
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

FAECAL VIRAL STUDIES (ADENOVIRUS ANTIGEN), Code: FVSA

Specimen: Stool (Sterile container with attached spatula)
 Collection: -
 Specimen Stability: 72 hours at 2-8°C
 Rejection Criteria: Specimens other than stool
 Method: Rapid strip test
 Useful for: Detection of Adenovirus in stool
 Reference range: Negative
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

FAECAL VIRAL STUDIES (ROTAVIRUS ANTIGEN), Code: FVSR

Specimen: Stool (Sterile container with attached spatula)
 Collection: -
 Specimen Stability: 72 hours at 2-8°C
 Rejection Criteria: Specimens other than stool
 Method: Rapid strip test
 Useful for: Detection of rotavirus in stool
 Reference range: Negative
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

FAECAL MICROSCOPY, CULTURE & SENSITIVITY Code: FMC

Specimen: Stool (Sterile container with attached spatula)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: Specimens other than stool
 Method: Conventional culture and light microscopy
 Useful for: Diagnosis of parasitic infection, *Salmonella*, *Shigella*, *Vibrio Cholerae* and *Campylobacter*.
 Reference range: Negative
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

FAECAL MICROSCOPY, Code: FM

Specimen: Stool (Sterile container with attached spatula)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: Specimens other than stool
 Method: Light microscopy
 Useful for: Detection of parasites in stool
 Reference range: No ova, cysts and parasitic infection
 Performed: Daily
 TAT: 1 day, STAT
 Section: Microbiology (ext 2436)

FERRITIN, Code: FER

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Screening test for iron status
 Reference range: See laboratory report
 Performed: Every Tuesday and Friday
 TAT: 3 days
 Section: Immunology (ext 2407)

FIBRINOGEN LEVEL, Code: FIB

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line indicated on tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, Haemolysed and/or clotted, lipaemic, icteric
 Method: Coagulometric Measurement
 Useful for: Aid in the diagnosis of fibrinogenaemia, disseminated, Intravascular coagulation and fibrinolysis
 Reference range: Refer Laboratory Report
 Performed: Daily, STAT
 TAT: 1 day
 Section: Haematology (ext 2408)

Refer to changes highlighted yellow

FLUID ALBUMIN, Code: ALBFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Colorimetric with Bromcresol Green
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID AMYLASE, Code: FLUAMY

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Enzymatic/Colorimetric with EPS
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID BIOCHEMISTRY – PROTEIN & GLUCOSE, Code: FLB

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

FLUID CHLORIDE, Code: CLFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID CREATININE, Code: CRTFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Kinetic Alkaline Picrate
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID CULTURE only, Code: FLC

Specimen: Body Fluid (Sterile screw-capped container, 1 ml)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: -
 Method: Conventional culture
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

FLUID FEME only (Gram Stain, Microscopy and cell count), Code: FLF

Specimen: Body Fluid (Sterile Screw-capped container, 1 ml)
 Collection: -
 Specimen Stability: Send to laboratory as soon as possible
 Rejection Criteria: -
 Method: Light microscopy
 Useful for: Presumptive diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 1 day, STAT
 Section: Microbiology (ext 2436)

FLUID LACTATE, Code: LACFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Enzymatic, Lactic acid to pyruvate
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID LACTATE DEHYDROGENASE (LDH), Code: LDHFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: IFCC, Lactate to Pyruvate
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID POTASSIUM, Code: POTFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID SODIUM , Code: NAFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FLUID TRIGLYCERIDE, Code: TGFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Glycerol Phosphate Oxidase
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

FLUID UREA, Code: URFLU

Specimen: Body fluid (Sterile container)
 Collection: Fresh sample and send to lab immediately
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Kinetic Urease
 Useful for: Identification of fluid effusions
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

FOLATE, Code: FOL

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Preferably fasting specimen. Send to the lab immediately.
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Investigation of megaloblastic anaemia and assessment of folate deficiency
 Reference range: See laboratory report
 Performed: Every Tuesday and Friday
 TAT: 3 days
 Section: Immunology (ext 2407)

FOLLICLE STIMULATING HORMONE, Code: FSH

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Evaluating the hypothalamic-pituitary-gonadal axis in diagnosing conditions such as amenorrhea, androgen deficiency, and gonadal dysfunction.**
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

FREE THYROXINE (FREE T4), Code: FT4

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of hyperthyroidism and hypothyroidism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Immunology (ext 2407)

FREE TRIIODOTHYRONINE (FREE T3), Code: FT3

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of hyperthyroidism and hypothyroidism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Immunology (ext 2407)

FULL BLOOD COUNT, Code: FBC

Specimen: Whole blood, (EDTA purple top)
 Collection: 4ml for EDTA tube and 500ul for microtainer
 Specimen Stability: 4 hours at room temperature
 Rejection Criteria: Haemolysed and/or clotted
 Method: Test includes machine operated differential count by light scattering flow cytometry
 Useful for: Provides information of cells in the blood
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

Refer to changes highlighted yellow

FUNGAL MICROSCOPY, Code: KOH

Specimen: Any (Sterile container)
 Collection: -
 Specimen Stability: -
 Rejection Criteria: -
 Method: Light microscopy
 Useful for: Detection of fungal infection
 Reference range: Negative
 Performed: Office hours only
 TAT: 24-48 hours
 Section: Microbiology (ext 2436)

FUNGAL MICROSCOPY AND CULTURE, Code: FUN

Specimen: Any (Sterile container)
 Collection: -
 Specimen Stability: -
 Rejection Criteria: -
 Method: Light microscopy & Conventional Culture
 Useful for: Detection of fungal infection
 Reference range: Negative
 Performed: Daily
 TAT: 4 weeks
 Section: Microbiology (ext 2436)

GAMMA-GLUTAMYL TRANSFERASE, Code: GGT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Enzymatic colorimetric (IFCC)
 Useful for: **Aids in diagnosis and treatment of liver diseases such as cirrhosis, primary and secondary liver tumours.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (Ext 2407)

Refer to changes highlighted yellow

GASTROINTESTINAL PCR, Code: GIP

Specimen: Stool
 Collection: In Cary Blair transport medium. Do not exceed the maximum filling line indicated on the label.
 Specimen Stability: 48 hours- Room temperature (20-25°C)/ 72 hours- Refrigeration (2-8°C)
 Rejection Criteria: Samples not in Cary Blair transport medium
 Method: BioFire FilmArray Torch
 Useful for: For the detection and identification of nucleic acids from multiple bacteria, viruses and parasites from individuals with signs and symptoms of gastrointestinal infection.
 Reference range: Not Detected
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

GENITAL SWAB CULTURE & SENSITIVITY, Code: GSC

Specimen: Genital Swab (Gel swab)
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab
 Method: Conventional culture and Gram Stain
 Useful for: Detection of yeast and bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

GENTAMICIN LEVEL (RANDOM), Code: GEN

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

GENTAMICIN LEVEL (PEAK), Code: GENP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Collect 30 min after end of IV infusion or 1 hour after IM injection
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

GENTAMICIN LEVEL (TROUGH), Code: GENT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Collect specimen immediately before next dose
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

GLOMERULAR FILTRATION RATE, Code: GFR

Tests include Creatinine and Estimated Glomerular Filtration Rate

Note: If serum creatinine was previously ordered, please use add-on test code: EGFR

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: **Please provide height (cm) for paediatric patients (< 18 years old)**
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Kinetic Alkaline Picrate
 Useful for: Renal function test
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (Ext 2407)

GLUCOSE 6-PHOSPHATE DEHYDROGENASE, Code: GPD

Specimen: Neonatal cord blood or whole blood (EDTA, purple top, 3mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Fluorescence Spot Test
 Useful for: Screening test for G6PD deficiency. (Note: any recent blood transfusion for the last 30 days or acute haemolysis can affect the results obtained with this test)
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (Ext 2407)

GLUCOSE FASTING, Code: GLF

Specimen: Plasma (Sodium fluoride, grey top, 3mL) preferred or serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Fasting for at least 8 hours
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: Not fasting or fasting for less than 8 hours.
 Method: Hexokinase/G-6-PDH
 Useful for: Diagnosis of diabetes mellitus
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (Ext 2407)

GLUCOSE 1 HR. POST PRANDIAL, Code: GLU1

Note: Please order this test together with GLUCOSE, FASTING (GLF)

Specimen: Plasma (Sodium fluoride, grey top, 3mL) preferred or serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: Sample collection time not labelled properly on tube
 Method: Hexokinase/G-6-PDH
 Useful for: Diagnosis of diabetes mellitus
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (Ext 2407)

Refer to changes highlighted yellow

GLUCOSE 2 HRS. POST PRANDIAL, Code: GLU2

Note: Please order this test together with GLUCOSE, FASTING (GLF)

Specimen: Plasma (Sodium fluoride, grey top, 3mL) preferred or serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)

Collection: NA

Specimen Stability: 7 days at 2-8°C

Rejection Criteria: Sample collection time not labelled properly on tube

Method: Hexokinase/G-6-PDH

Useful for: Diagnosis of diabetes mellitus

Reference range: See laboratory report

Performed: Daily

TAT: 1 day

Section: Biochemistry (Ext 2407)

GLUCOSE RANDOM, Code: GLR

Specimen: Plasma (Sodium fluoride, grey top, 3mL) preferred or serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)

Collection: NA

Specimen Stability: 7 days at 2-8°C

Rejection Criteria: NA

Method: Hexokinase/G-6-PDH

Useful for: Diagnosis of diabetes mellitus

Reference range: See laboratory report

Performed: Daily

TAT: 1 day, STAT

Section: Biochemistry (Ext 2407)

GLUCOSE TOLERANCE TEST (2 SPECIMENS), Code: GTT2

Specimen: Plasma (Sodium fluoride, grey top, 3mL) and Urine (Sterile container).

Collection: Submit 2 groups of specimens:
 1. Fasting (plasma and urine)
 2. 2 hours after glucose (75g) intake (plasma and urine)

Specimen Stability: 7 days at 2-8°C

Rejection Criteria: Sample collection time not labelled properly (fasting/ 2nd hour) on tubes

Method: Hexokinase/G-6-PDH (plasma), Glucose oxidase/peroxidase (urine)

Useful for: Diagnosis of diabetes mellitus

Reference range: See laboratory report

Performed: Daily

TAT: 1 day

Section: Biochemistry (Ext 2407)

GLYCATED HAEMOGLOBIN (HbA1C), Code: HBA1C

Specimen: Whole blood (EDTA, purple top, 3mL).
 Collection: If BG and/or FBC is also requested, please take an extra EDTA purple top for this test.
 Mix specimen well to avoid clotting
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: Clotted specimen
 Method: Enzymatic
 Useful for: Long term monitoring of glucose control in diabetes mellitus
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 1 day
 Section: Biochemistry (ext 2407)

GRAM STAIN MICROSCOPY, Code: GRAM

Specimen: Any (Sterile container or Swab)
 Collection: -
 Specimen Stability: -
 Rejection Criteria: -
 Method: Conventional Gram Stain
 Useful for: -
 Reference range: -
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

HDL CHOLESTEROL, Code: HDL

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: N/A
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: Grossly haemolysed, Icteric or Lipemic
 Method: Colorimetric with accelerator selective detergent
 Useful for: Cardiovascular risk assessment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

HEPATITIS A IGM ANTIBODY, Code: HAM

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of acute or recent hepatitis A infection
 Reference range: Non-reactive
 Performed: Office hours only
 TAT: 1-7 days
 Section: Virology and Serology (ext 2322)

HEPATITIS A IGG ANTIBODY, Code: HAG

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of Hepatitis A infection or vaccination
 Presence of IgG anti-HAV, with a non-reactive IgM anti-HAV test result, implies past infection with hepatitis A virus (HAV) or vaccination against HAV
 Reference range: Non-reactive or Reactive (if immunized)
 Performed: Office hours only
 TAT: 1-7 days
 Section: Virology and Serology (ext 2322)

HEPATITIS B CORE TOTAL ANTIBODY, Code: HBT

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Detection for acute or chronic Hepatitis B infection, persists for life
 Reference range: Non-reactive
 Performed: Office hours only
 TAT: 1-3 days
 Section: Virology and Serology (ext 2322)

Refer to changes highlighted yellow

HEPATITIS B DNA VIRAL LOAD, Code: HBVDNA

Specimen: Plasma (EDTA, purple top, 4 mL)
 Collection: NA
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: NA
 Method: NAAT
 Useful for: Detection and quantification of hepatitis B virus (HBV) DNA in patients with HBV infection (ie, HBsAg-positive). Monitor response to Anti-HBV therapy.
 Reference range: Not Detected
 Performed: Office Hours Only
 TAT: 1-3 days
 Section: Virology and Serology (ext 2322)

HEPATITIS B SURFACE ANTIGEN, Code: HB1

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of acute or chronic Hepatitis B infection; earliest routine indicator.
Repeatedly reactive result will be automatically reflexed to HBsAg Confirmatory (Neutralisation) assay.
 Reference range: Non-reactive
 Performed: Office hours only
 TAT: 1-2 days
 Section: Virology and Serology (ext 2322)

HEPATITIS B SURFACE ANTIBODY, Code: HB2

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Used to detect natural immunity or vaccination to Hepatitis B
 Reference range: ≥ 10 mIU/mL indicates evidence of protective immunity either from past infection or vaccination
 Performed: Office hours only
 TAT: 1 day
 Section: Virology and Serology (ext 2322)

Refer to changes highlighted yellow

HEPATITIS C ANTIBODY, Code: HCV

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Detects acute or chronic infection to Hepatitis C
 Repeatedly reactive result should be confirmed by Hepatitis C confirmatory assay.
 Reference range: Non-reactive
 Performed: Office hours only
 TAT: 1 day
 Section: Virology and Serology (ext 2322)

HIGH-SENSITIVE C-REACTIVE PROTEIN, Code: HSCRP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 2 months at 2-8°C
 Rejection Criteria: **NA**
 Method: Turbidimetric/Immunturbidimetric
 Useful for: **A diagnostic indicator of infections and inflammation particularly in paediatric patients and monitoring response to therapy. Useful for risk management of coronary heart disease.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Biochemistry (ext 2407)

HIV 1&2 AG/AB, Code: HIV

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Determination of HIV p24 antigen and antibodies to HIV -1 and HIV - 2
 Repeatedly reactive result should be confirmed by HIV confirmatory assay.
 Reference range: Non-reactive
 Performed: Office hours only
 TAT: 1-2 days
 Section: Virology and Serology (ext 2322)

Refer to changes highlighted yellow

HUMAN CHORIONIC GONADOTROPIN BETA TOTAL (β -hCG), Code: QUA

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Early detection and monitoring of pregnancy. A tumour marker for hydatiform mole, choriocarcinoma and testicular cancer
 Reference range: See laboratory report
 Performed: Office hours only
 TAT: 2 days, STAT
 Section: Immunology (ext 2407)

INFERTILITY STUDIES (FEMALE), Code: ISF

Tests include Estradiol, FSH, LH, Progesterone, Prolactin and Beta-HCG

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Investigation of infertility for females
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

INFERTILITY STUDIES (MALE), Code: IFM

Tests include FSH, LH, Prolactin and Testosterone

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Investigation of infertility for males
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

INSULIN , Code: INS	
Specimen:	Serum (SST, yellow top, 5mL) or Plasma (EDTA, purple top, 4mL)
Collection:	Fasting for at least 8 hours. Send the sample to the lab immediately.
Specimen Stability:	7 days at -10°C or colder
Rejection Criteria:	Not fasting or fasting for less than 8 hours.
Method:	Chemiluminescent Microparticle Immunoassay (CMIA)
Useful for:	Aids in diagnosis of diabetes mellitus, differential diagnosis of fasting hypoglycaemia (factitious hypoglycaemia, insulin autoimmune hypoglycaemia and insulinoma)
Reference range:	See laboratory report
Performed:	Monday & Thursday
TAT:	2-3 days
Section:	Immunology (ext 2407)
IRON, Code: FE	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
Collection:	NA
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	Colorimetric with Ferene
Useful for:	Evaluation of Iron metabolism. Aid in diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis and chronic renal disease.
Reference range:	See laboratory report
Performed:	Every Tuesday and Friday
TAT:	3 days
Section:	Biochemistry (ext 2407)
IRON STUDIES, Code: IS	
<i>Tests include Total Iron, Ferritin, Transferrin, Transferrin Saturation and TIBC</i>	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
Collection:	NA
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	Refer to individual test method
Useful for:	Evaluation of Iron metabolism. Aid in diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis and chronic renal disease.
Reference range:	See laboratory report
Performed:	Every Tuesday and Friday
TAT:	3 days
Section:	Biochemistry (ext 2407)

Refer to changes highlighted yellow

INTERNATIONAL NORMALISED RATIO, Code: INR

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line indicated on tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Derived from calculation
 Useful for: Monitoring of warfarin dosage
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

LACTATE DEHYDROGENASE, Code: LDH

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: NA
 Method: IFCC, lactate to pyruvate
 Useful for: Non-specific marker of cellular damage. Aid in the differential diagnosis and treatment of hemolytic anaemia, liver diseases, cardiac diseases, and used as a tumour marker in some malignancies.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

LACTIC ACID/ LACTATE, Code: LAC

Specimen: Plasma (Sodium fluoride, grey top, 3 mL)
 Collection: Avoid the use of tourniquet if possible and send to the Lab immediately within 30 minutes for rapid centrifugation
 Specimen Stability: 3 days at room temperature
 Rejection Criteria: NA
 Method: Enzymatic, Lactic acid to pyruvate
 Useful for: Evaluation of metabolic and lactic acidosis
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

LIPID PROFILE (CORONARY RISK FACTORS), Code: LIP

Tests include Total Cholesterol, LDL Cholesterol, HDL Cholesterol, Triglycerides and Total Cholesterol/HDL ratio

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: **Fasting is recommended for at least 8 hours**
 Specimen Stability: 5 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Refer to individual test method
 Useful for: Diagnosis of hyperlipidemia and cardiovascular risk assessment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

LIVER FUNCTION TEST, Code: LFT

Tests include Total Protein, Albumin, Globulin, A/G Ratio, Bilirubin Total, ALP, AST/SGOT, ALT/SGPT and GGT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Refer to individual test method
 Useful for: **Assessment of liver function**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

LUTEINIZING HORMONE, Code: LH

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Evaluating the hypothalamic-pituitary-gonadal axis in diagnosing conditions such as amenorrhea, androgen deficiency, and gonadal dysfunction.**
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

MAGNESIUM, Code: MG

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Enzymatic
 Useful for: Diagnosis and monitoring of hypo- and hypermagnesemia
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

MAGNESIUM (24 HOURS URINE), Code: UMG24

Specimen: Urine (24hr urine container)
 Collection: Preserve urine using 20 to 30 mL of 6M HCL and store refrigerated during collection period
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: NA
 Method: Enzymatic
 Useful for: Diagnosis and monitoring of hypo- and hypermagnesemia
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

MALARIAL PARASITES, Code: MS

Specimen: Whole blood, (EDTA purple top, 4ml)
 Collection: Preferably at peak of fever
 Specimen Stability: 48 hours from the time of collection
 Rejection Criteria: Haemolysed and/or clotted
 Method: Field stain, Microscopic analysis
 Useful for: Detection and identification of malarial parasites
 Reference range: No malaria parasite detected
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

Refer to changes highlighted yellow

MENINGITIS/ ENCEPHALITIS PCR, Code: MEP

Specimen: Cerebro-spinal Fluid (CSF)
 Collection: Minimum volume – 0.2 mL (200 µL)
 Specimen Stability: 24 hours – Room temperature (Approx. 23°C)/ 7 days- Refrigeration (Approx. 4°C)
 Rejection Criteria: -
 Method: BioFire FilmArray Torch
 Useful for: For the detection and identification of bacterial, viral, and yeasts nucleic acids from multiple meningitis and encephalitis pathogens
 Reference range: Not Detected
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

METHOTREXATE, Code: MTX

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Protect sample from light and send to lab immediately
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **Sample not protected from light**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Immunology (ext 2407)

MISCELLANEOUS CULTURE & SENSITIVITY, Code: RSCM

Specimen: Miscellaneous Swab (Gel swab), Tissue (Sterile Container)
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab, Tissue in formalin, Non-sterile container
 Method: Conventional culture and Gram Stain
 Useful for: To detect bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

MICROALBUMIN (URINE RANDOM), Code: MAS

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Turbidimetric/Immunoturbidimetric
 Useful for: Early detection of diabetic nephropathy
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

MICROALBUMIN (24 HOURS URINE), Code: UALB24

Specimen: Urine (24hr urine container)
 Collection: Preserve urine using 20 to 30 mL of 6M HCL and store refrigerated during collection period
 Specimen Stability: Assay immediately
 Rejection Criteria: NA
 Method: Turbidimetric/Immunoturbidimetric
 Useful for: Early detection of diabetic nephropathy
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

MICROFILARIA PARASITES, Code: MFS

Specimen: Whole blood, (EDTA purple top, 4ml)
 Collection: Send immediately to the laboratory
 Specimen Stability: 48 hours from the time of collection
 Rejection Criteria: Haemolysed and/or clotted
 Method: Field stain, Microscopic analysis
 Useful for: Detection of microfilariae in blood to aid in the diagnosis of Filariasis
 Reference range: Negative
 Performed: Office hours only
 TAT: 1 day
 Section: Haematology (ext 2408)

Refer to changes highlighted yellow

MRSA SWAB CULTURE & SENSITIVITY, Code: Refer note below

Specimen: Any (Gel swab), indicate source of specimen
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab, Source of specimen not indicated if more than one swab
 Method: Conventional culture
 Useful for: To detect Meticillin-resistant *Staphylococcus aureus*
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

NOTE: Code for MRSA Swab culture & sensitivity is MSA(no. of swabs)
 (i.e. For one MRSA swab, code is MSA1. For two MRSA swabs, code is MSA2....)

M. TUBERCULOSIS COMPLEX / RIFAMPICIN RESISTANCE PCR, Code: MTBR

Specimen: Sputum (3 mL), or CSF (2 mL)
 Collection: Collect in a dedicated sterile screw-capped container, and deliver to the lab immediately
 Specimen Stability: Room temperature, 72 hours
 Rejection Criteria: -
 Method: Xpert MTB/RIF Assay using Real-time Polymerase Chain Reaction
 Useful for: Rapid detection of *Mycobacterium tuberculosis* DNA for the diagnosis of pulmonary tuberculosis, and presumptive detection of rifampin resistance based on the presence of resistance-associated mutations
(This test should always be performed in conjunction with mycobacterial culture)
 Reference range: Not Detected
 Performed: Office hours only
 TAT: 1-3 days
 Section: Virology and Serology (ext 2322)

Refer to changes highlighted yellow

NEONATAL BILIRUBIN Code: NBI

Tests include Bilirubin Total & Bilirubin Direct (Conjugated)

Specimen: Serum (Microtainer SST, 500uL) or plasma (Microtainer Li-Heparin, 500uL)
 Collection: Send to the Lab immediately
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Diazonium salt
 Useful for: **Diagnosing neonatal jaundice and assessing the effectiveness of phototherapy in neonates.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day, STAT**
 Section: Biochemistry (ext 2407)

NOSE SWAB (NASAL CULTURE & SENSITIVITY), Code: RSCN

Specimen: Nose Swab (Gel swab)
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab
 Method: Conventional culture and Gram Stain
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

N-TERMINAL PRO B-TYPE NATRIURETIC PEPTIDE, Code: NTPROBNP

Specimen: Serum (SST, yellow top, 5mL) **or plasma (Li-Heparin green top, 4 mL) or Plasma (EDTA, purple top, 4mL)**
 Collection: NA
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Aids in diagnosis of congestive heart failure (CHF), detection of mild forms of cardiac dysfunction and assessment of heart failure severity. Useful for monitoring treatment in patients with left ventricular dysfunction.
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

OSMOLALITY (SERUM), Code: OSS

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Freezing point osmometry
 Useful for: Assessment of fluid and electrolyte balance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

OSMOLALITY (URINE), Code: UOS

Specimen: Urine (Sterile container)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Freezing point osmometry
 Useful for: Assessment of fluid and electrolyte balance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

PARATHYROID HORMONE , Code: PTH

Specimen: **Plasma (EDTA, purple top, 4mL) or plasma (Li-Heparin green top, 4 mL)**
 Collection: NA
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Aids in differential diagnosis of hypercalcemia, hypocalcemia, and parathyroid disorders. Useful in monitoring dialysis patients to manage renal osteodystrophy.
 Reference range: See laboratory report
 Performed: Monday, Wednesday & Friday
 TAT: 2 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

PHENOBARBITAL, Code: PHNO3

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 6 months at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: **Therapeutic Drug Monitoring**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

PHENYTOIN, Code: PHNY2

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 1 month at 2-8°C
 Rejection Criteria: **NA**
 Method: Enzyme Immunoassay
 Useful for: **Therapeutic Drug Monitoring**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

PHOSPHATE/PHOSPHORUS, Code: PO4

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: **3 days at 2-8°C**
 Rejection Criteria: **NA**
 Method: Phosphomolybdate Formation
 Useful for: **Used in the diagnosis and treatment of various disorders including parathyroid gland, kidney diseases and vitamin D imbalance.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

PHOSPHATE/PHOSPHORUS (URINE RANDOM), Code: UPO4

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Phosphomolybdate Formation
 Useful for: Used in the diagnosis and treatment of various disorders including parathyroid gland, kidney diseases and vitamin D imbalance.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

PHOSPHATE/PHOSPHORUS (24 HOURS URINE), Code: UPO424

Specimen: Urine (24hr urine container)
 Collection: Preserve urine using 20 to 30 mL of 6M HCL and store refrigerated during collection period
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Phosphomolybdate Formation
 Useful for: Used in the diagnosis and treatment of various disorders including parathyroid gland, kidney diseases and vitamin D imbalance.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

POTASSIUM, Code: POT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

POTASSIUM (URINE RANDOM), Code: UPOT

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 2 months at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation / assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

POTASSIUM (24 HOURS URINE), Code: UPOT24

Specimen: Urine (24hr urine container)
 Collection: No preservative needed and store refrigerated during collection period
 Specimen Stability: 2 months at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation / assessment of electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

PROCALCITONIN, Code: PCT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Aid in differentiation of bacterial infection, assessment of severity of bacterial infection in suspect sepsis patients, aid in prognosis and monitoring antibiotic therapy.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: **1 day**
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

PROGESTERONE, Code: PRG

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 10 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Evaluation of ovarian function and abnormal pregnancy. Detection of progesterone-secreting tumour. Aid in optimizing IVF procedure for successful embryo implantation.**
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days, STAT
 Section: Immunology (ext 2407)

PROLACTIN, Code: PRL

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Evaluation of subfertility, hypogonadism and pituitary gland function
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

PROSTATE-SPECIFIC ANTIGEN TOTAL, Code: PSA

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 24 hours at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Tumour marker for prostate cancer. Indicator for benign hypertrophy and prostatitis.**
 Reference range: See laboratory report
 Performed: **Office hours only**
 TAT: **2 days**
 Section: Immunology (ext 2407)

PROTEIN (URINE RANDOM), Code: UTP

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Turbidimetric with Benzethonium Chloride
 Useful for: Indicator of renal impairment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

PROTEIN (24 HOURS URINE), Code: UTP24

Specimen: Urine (24hr urine container)
 Collection: No preservative needed and store refrigerated during collection period
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: NA
 Method: Turbidimetric with Benzethonium Chloride
 Useful for: Indicator of renal impairment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

PROTEIN CREATININE RATIO, Code: PCRTR

Specimen: Urine (Sterile container or 24hr urine container)
 Collection: Random or 24 hours. No preservative needed
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Indicator of renal impairment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

PROTEIN TOTAL, Code: TP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at room temperature and at 2-8°C
 Rejection Criteria: **NA**
 Method: Colorimetric with Biuret
 Useful for: **Aid in diagnosis and treatment of various diseases involving liver, kidney or bone marrow, other metabolic or nutritional disorders.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

PROTHROMBIN TIME, Code: PT

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line as indicated on tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Coagulometric Measurement
 Useful for: Screening test for clotting disorders. Monitoring of anticoagulation therapy
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

PT 50% CORRECTION, Code: PT50

Specimen: Plasma (Sodium Citrate, blue top)
 Collection: Up to line as indicated on tube
 Specimen Stability: 2 hours from the time of collection
 Rejection Criteria: Under filled, over filled, haemolysed and/or clotted, lipaemic, icteric
 Method: Coagulometric Measurement
 Useful for: Follow up test to investigate cause of unexplained prolonged PT. Differentiates between possibility of factor deficiency or acquired coagulation inhibitors of the extrinsic and common pathways.
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day
 Section: Haematology (ext 2408)

Refer to changes highlighted yellow

PUS SWAB CULTURE & SENSITIVITY, Code: PUS

Specimen: Pus (Gel Swab/Sterile container/Syringe with needle removed)
 Collection: Indicate source of specimen
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry swab, samples with formalin
 Method: Conventional culture and Gram Stain
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

RAPID PLASMA REAGIN, Code: RPR

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 5 days at 2-8°C
 Rejection Criteria: NA
 Method: Flocculation test
 Useful for: Detection of reagin antibodies associated with syphilis; to monitor response to syphilis treatment
 Reference range: NA
 Performed: Office hours only
 TAT: 2 days
 Section: Virology and Serology (ext 2322)

RENAL FUNCTION TEST, Code: RFT

Tests include Electrolytes, Urea, Creatinine and Uric Acid

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: NA
 Method: See individual test method
 Useful for: Assessment of renal function
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

RESPIRATORY PCR PANEL, Code: RP	
Specimen:	Nasopharyngeal swab (NPS)
Collection:	Immediately place swab in the medium provided Minimum volume – 0.3 mL (300 µL)
Specimen Stability:	4 hours – Room temperature (15-25°C)/ 3 days- Refrigeration (2-8°C)/ 30 days – Frozen (≤-15 °C or ≤-70°C)
Rejection Criteria:	-
Method:	BioFire FilmArray Torch
Useful for:	For qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) obtained from individuals suspected of respiratory tract infections
Reference range:	Not Detected
Performed:	Daily
TAT:	1 day
Section:	Microbiology (ext 2436)
RESPIRATORY SPECIMEN CULTURE & SENSITIVITY, Code: RSCS	
Specimen:	Sputum, Bronchoalveolar Lavage (BAL), Endotracheal Aspirate (ETA), Nasopharyngeal Aspirate (NPA)
Collection:	-
Specimen Stability:	Room temperature, 24 hours
Rejection Criteria:	Sputum visually salivary
Method:	Conventional culture and Gram Stain
Useful for:	Diagnosis of bacterial infection
Reference range:	-
Performed:	Daily
TAT:	48-72 hours
Section:	Microbiology (ext 2436)
RESPIRATORY VIRUS (FLU A/B, RSV, RESP. ADENOVIRUS) SCREENING, Code: RES3	
Specimen:	Nasal/Nasopharyngeal swab (Dry Swab)
Collection:	NA
Specimen Stability:	8 hours at 2-8°C
Rejection Criteria:	Gel swab
Method:	Rapid immunochromatographic assay
Useful for:	Determination of type of Influenza infection (type A or B), Respiratory Syncytial virus and Adenovirus antigens
Reference range:	Negative
Performed:	Daily
TAT:	1 day, STAT
Section:	Virology and Serology (ext 2322)

RETICULOCYTE COUNT, Code: RET

Specimen: Whole blood (EDTA purple top)
 Collection: 4ml for EDTA tube and 500ul for microtainer
 Specimen Stability: 4 hours at room temperature
 Rejection Criteria: Haemolysed and/or clotted
 Method: Supravital staining, manual technique or by electronic blood cell analyser
 Useful for: Assessment of erythropoietic activity and evaluation of anemia
 Reference range: See Laboratory Report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Haematology (ext 2408)

RHEUMATOID FACTOR, Code: RHE

Specimen: Serum (SST yellow top, 5mL)
 Collection: NA
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Turbidimetric/Immuno-turbidimetric
 Useful for: Supports diagnosis of Rheumatoid Arthritis and evaluation of ankylosing spondylitis, sjogren's syndrome, scleroderma, dermatomyositis and SLE
 Reference range: See laboratory report
 Performed: Every Monday, Wednesday and Friday
 TAT: 2 days
 Section: Biochemistry (ext 2407)

RICKETTSIAL SEROLOGY (Weil Felix Screen), Code: RIC

Specimen: Serum (SST tube, Yellow top, 1 mL)
 Collection: -
 Specimen Stability: 2°C to 8°C for 48 hours
 Rejection Criteria: Grossly Haemolysed/ Grossly Lipaemic/Contaminated
 Method: Agglutination
 Useful for: Diagnosis of rickettsial infection
 Reference range: Titre < 1:80
 Performed: Office hours only
 TAT: 1 day
 Section: Microbiology (ext 2436)

RISK STRATIFICATION USING TROPONIN-I, Code: RISK

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 24 hours at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: To aid in stratifying the risk of cardiovascular disease, including cardiovascular death, MI, coronary revascularization, heart failure, or ischemic stroke in asymptomatic individuals
 Reference range: See laboratory report
 Performed: **Daily**
 TAT: 1 day
 Section: Biochemistry (ext 2407)

RUBELLA IGG, Code: RUG

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 14 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Detection of antibodies and immunity to rubella
 Reference range: ≥ 10 mIU/mL indicates evidence of protective immunity either from past infection or vaccination
 Performed: Office hours only
 TAT: **7-14** days
 Section: Virology and Serology (ext 2322)

SALMONELLA TYPHI IgG/IgM, CODE: STGM

Specimen: Serum or Plasma **or Whole Blood**
 Collection: **SST or Plain Tube/ Lithium Heparin Tube/ K2-EDTA Tube**
 Specimen Stability: 72 hours at 2-8°C (Serum/ Plasma)/ **48 hours at 2-8°C (Whole Blood)**
 Rejection Criteria: **Haemolysed**
 Method: Immunochromatographic Test
 Useful for: Detection of Salmonella typhi IgG and IgM antibodies
 Reference range: Negative
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

SEMEN CULTURE & SENSITIVITY, Code: SEMEN

Specimen: Fresh Semen (Sterile container)
 Collection: -
 Specimen Stability: Send to lab as soon as possible
 Rejection Criteria: -
 Method: Conventional culture and Gram Stain
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

SODIUM, Code: NA

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 2 weeks at 2-8°C
 Rejection Criteria: **NA**
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation of fluid and electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

SODIUM (URINE RANDOM), Code: UNA

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 45 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation of fluid and electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

SODIUM (24 HOURS URINE), Code: UNA24

Specimen: Urine (24hr urine container)
 Collection: No preservative needed and store refrigerated during collection period
 Specimen Stability: 45 days at 2-8°C
 Rejection Criteria: NA
 Method: Indirect Ion-selective electrode potentiometry
 Useful for: Evaluation of fluid and electrolyte imbalance
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

STERILITY TEST (ENVIRONMENTAL), Code: *Refer to note below*

Specimen: Any (Agar plate)
 Collection: -
 Specimen Stability: -
 Rejection Criteria: -
 Method: Open plate technique/conventional culture
 Useful for: Detection of pathogens
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

NOTE: Code for sterility test is STER(no. of specimen)
(i.e. For one specimen, code is STER1. For two specimen, code is STER2. And so on and so forth).

SYNOVIAL FLUID for CRYSTAL EXAMINATION, Code: SFC

Specimen: Synovial Fluid (Sterile container)
 Collection: -
 Specimen Stability: -
 Rejection Criteria: -
 Method: Wet preparation for microscopic analysis
 Useful for: Detection of crystals and specific infections
 Reference range: -
 Performed: Office hours only
 TAT: 1 day
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

SYPHILIS SCREENING, Code: SYP

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Determination of anti-Treponema pallidum antibodies in human Serum (Recommended for syphilis screening)
 Reference range: Non-reactive
 Performed: Office hours only
 TAT: 1-2 days
 Section: Virology and Serology (ext 2322)

TACROLIMUS, Code: TACRO

Specimen: Whole blood (EDTA, purple top, 3mL).
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Therapeutic Drug Monitoring for Organ Transplant
 Reference range: See laboratory report
 Performed: Daily, STAT
 TAT: 1 day
 Section: Immunology (ext 2407)

TESTOSTERONE TOTAL, Code: TES

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Evaluation of subfertility in males; hirsutism and virilisation in females
 Reference range: See laboratory report
 Performed: Every Monday and Thursday
 TAT: 3 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

THEOPHYLLINE, Code: THE

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 3 months at 2-8°C
 Rejection Criteria: **NA**
 Method: Enzyme Immunoassay
 Useful for: Monitoring of theophylline
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

THYROGLOBULIN , Code: THY

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL) or Plasma (EDTA, purple top, 4mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Used as tumour marker for papillary and follicular thyroid cancer. Aids in diagnosis of Thyrotoxicosis factitial and monitoring after total thyroid ablation (thyroidectomy)
 Reference range: See laboratory report
 Performed: Wednesday & Saturday
 TAT: 2-3 days
 Section: Immunology (ext 2407)

THYROID AUTOANTIBODIES, Code: THA

Tests include Anti-Thyroid Receptor Antibody, Anti-Thyroid Peroxidase, Anti-Thyroglobulin

Specimen: Serum (SST, yellow top, 5mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: NA
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Aids in differential diagnosis of autoimmune thyroid diseases.
 Reference range: See laboratory report
 Performed: Wednesday & Saturday
 TAT: 2-3 days
 Section: Immunology (ext 2407)

Refer to changes highlighted yellow

THYROID FUNCTION TEST, Code: TFT

Tests include FT4, FT3 and Thyroid Stimulating Hormone

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 6 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of hyperthyroidism and hypothyroidism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Immunology (ext 2407)

THYROID STIMULATING HORMONE, Code: TSH

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Diagnosis of hyperthyroidism and hypothyroidism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Immunology (ext 2407)

THROAT SWAB (Culture and Sensitivity), Code: RSCT

Specimen: Throat Swab (Gel swab)
 Collection: -
 Specimen Stability: Room temperature, 24 hours
 Rejection Criteria: Dry Swab
 Method: Conventional culture
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

TIP CULTURE & SENSITIVITY, Code: RSCTIP

Specimen: IUCD, Catheter Tip, ETC, PICC, CVC, any tips (Sterile container)
 Collection: Indicate source of specimen
 Specimen Stability: 24 hours at room temperature
 Rejection Criteria: -
 Method: Conventional culture
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

TRANSFERRIN, Code: TRANS

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 3 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Immunospectrometric
 Useful for: Screening test for iron status
 Reference range: See laboratory report
 Performed: Every Tuesday and Friday
 TAT: 3 days
 Section: Biochemistry (ext 2407)

TREPONEMA PALLIDUM PARTICLE AGGLUTINATION, Code: TPA

Specimen: Serum (SST, Yellow top, 5 mL)
 Collection: NA
 Specimen Stability: 2 weeks at 2-8°C
 Rejection Criteria: Haemolysis, lipaemic, contaminated
 Method: Passive Particle Agglutination
 Useful for: Detection and titration of antibodies against causative agent of Syphilis (Recommended for syphilis confirmatory test)
 Reference range: Negative
 Performed: Office hours only
 TAT: 3 days
 Section: Virology and Serology (ext 2322)

TRIGLYCERIDE, Code: TG

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Fasting is recommended for 8-12 hours
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Enzymatic with Glycerol Phosphate Oxidase
 Useful for: Diagnosis of hyperlipidemia and cardiovascular risk assessment
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

TROPONIN-I HIGH SENSITIVE STAT, Code: TRO

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 24 hours at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: **Aid in diagnosis of myocardial infarction and assessment of 30-day and 90-day prognosis to all-cause mortality and major adverse cardiac events.**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Immunology (ext 2407)

TYPHOID ANTIBODY (Widal test), Code: TYP

Specimen: Serum (SST tube, Yellow top, 1 mL)
 Collection: -
 Specimen Stability: 2°C to 8°C for 48 hours
 Rejection Criteria: Grossly Haemolysed/ Grossly Lipaemic/Contaminated
 Method: Agglutination
 Useful for: Diagnosis of *Salmonella* infection
 Reference range: Titre < 1:80
 Performed: Office hours only
 TAT: 1 day
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

UREA, Code: UR	
Specimen:	Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
Collection:	NA
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	Kinetic
Useful for:	Assessment of renal function and differential diagnosis of prerenal, renal and postrenal hyperuremia.
Reference range:	See laboratory report
Performed:	Daily
TAT:	1 day, STAT
Section:	Biochemistry (ext 2407)
UREA (URINE RANDOM), Code: URU	
Specimen:	Urine (Sterile container)
Collection:	No preservative needed
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	Kinetic
Useful for:	Assessment of fluid balance and renal function
Reference range:	See laboratory report
Performed:	Daily
TAT:	1 day
Section:	Biochemistry (ext 2407)
UREA (24 HOURS URINE), Code: URU24	
Specimen:	Urine (24hr urine container)
Collection:	No preservative needed
Specimen Stability:	7 days at 2-8°C
Rejection Criteria:	NA
Method:	Kinetic
Useful for:	Assessment of fluid balance and renal function
Reference range:	See laboratory report
Performed:	Daily
TAT:	1 day
Section:	Biochemistry (ext 2407)

Refer to changes highlighted yellow

URIC ACID, Code: SUA

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 8 hours at room temperature, 3 days at 2-8°C
 Rejection Criteria: NA
 Method: Enzymatic colorimetric with Uricase
 Useful for: Aid in diagnosis and treatment of renal and metabolic disorders, including renal failure, gout, leukemia, psoriasis, diabetes, hypothyroidism, and atherosclerosis.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

URIC ACID (URINE RANDOM), Code: UUA

Specimen: Urine (Sterile container)
 Collection: No preservative needed
 Specimen Stability: 2 days at room temperature and at 2-8°C
 Rejection Criteria: NA
 Method: Enzymatic colorimetric with Uricase
 Useful for: Evaluation of uric acid metabolism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

URIC ACID (24 HOURS URINE), Code: UUA24

Specimen: Urine (24hr urine container)
 Collection: No preservative needed
 Specimen Stability: 2 days at room temperature and at 2-8°C
 Rejection Criteria: NA
 Method: Enzymatic colorimetric with Uricase
 Useful for: Evaluation of uric acid metabolism
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

URIC ACID/CREATININE RATIO, Code: UUACRTR

Specimen: Urine (Sterile container or 24hr urine container)
 Collection: No preservative needed
 Specimen Stability: 2 days at room temperature and at 2-8°C
 Rejection Criteria: NA
 Method: Refer to individual test method
 Useful for: Aid in diagnosis and treatment of renal and metabolic disorders, including renal failure, gout, leukemia, psoriasis, diabetes, hypothyroidism, and atherosclerosis.
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

URINE DRUG SCREEN, Code: UDRUG

(Tests include Amphetamine/Methamphetamine, Benzodiazepine, Cannabinoids, Cocaine, Ecstasy (MDMA), Opiates, Creatinine, Specific Gravity)

Specimen: Urine (Urine drug test container with temperature strip)
 Collection: Collect 2 specimens with at least 20 mL each. Seal the containers properly, date, time and initial on the seal stickers.
 Specimen Stability: 5 days at 2-8°C
 Rejection Criteria: Only one urine specimen received, not sealed, tampered container seal, volume is less than 20 mL
 Method: Enzyme Immunoassay
 Useful for: Qualitative detection of drug abuse
 Reference range: See laboratory report
 Performed: Every Tuesday, Thursday, and Saturday
 TAT: 2 days
 Section: Biochemistry (ext 2407)

URINE (Bilharzia/Helminths), Code: UBH

Specimen: Random Urine (Sterile container)
 Collection: -
 Specimen Stability: Up to 24 hours at 2-8°C
 Rejection Criteria: Specimen exceeds stability, container with preservative, non-sterile container
 Method: Wet slide preparation & Light Microscopy
 Useful for: Detection of Bilharzia or helminths in urine
 Reference range: -
 Performed: Office hours only
 TAT: 1 day
 Section: Microbiology (ext 2436)

Refer to changes highlighted yellow

URINE FEME, Code: MU

Specimen: Midstream Urine (Sterile container)
 Collection: -
 Specimen Stability: Up to 4 hours at room temperature; up to 24 hours at 2-8°C
 Rejection Criteria: Specimen exceeds stability, container with preservative, non-sterile container
 Method: iRICELL 2000
 Useful for: Urine analysis and microscopy
 Reference range: -
 Performed: Daily
 TAT: 1 day, STAT
 Section: Microbiology (ext 2436)

URINE FEME, CULTURE & SENSITIVITY Code: ZUCS

Specimen: Midstream Urine (Sterile container)
 Collection: -
 Specimen Stability: Up to 4 hours at room temperature; up to 24 hours at 2-8°C
 Rejection Criteria: Specimen exceeds stability, container with preservative, non-sterile container
 Method: Conventional culture & iRICELL 2000
 Useful for: Detection of urinary tract infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

URINE GLUCOSE, Code: UGLU

Specimen: Random Urine (Sterile container)
 Collection: A mid-stream urine specimen is recommended. **Avoid collecting during menstrual period as it may contaminate urine**
 Specimen Stability: Up to 2 hours at room temperature, up to 24 hours at 2-8°C
 Rejection Criteria: **Specimen exceeds stability**
 Method: Urine Dipstick (Glucose Oxidase)
 Useful for: Screening or management of diabetes mellitus
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

URINE KETONE, Code: UK

Specimen: Random Urine (Sterile container)
 Collection: A mid-stream urine specimen is recommended. **Avoid collecting during menstrual period as it may contaminate urine**
 Specimen Stability: Up to 2 hours at room temperature, up to 24 hours at 2-8°C
 Rejection Criteria: **Specimen exceeds stability**
 Method: Urine Dipstick (Sodium nitroprusside)
 Useful for: Screening for presence of ketoacidosis
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day
 Section: Biochemistry (ext 2407)

URINE RED CELL MORPHOLOGY, Code: RCM

Specimen: Random Urine (Sterile container)
 Collection: -
 Specimen Stability: Up to 2 hours at room temperature
 Rejection Criteria: Specimen exceeds stability, container with preservative
 Method: Phase contrast microscopy
 Useful for: Evaluation of Glomerular diseases
 Reference range: -
 Performed: Daily
 TAT: 1 day
 Section: Microbiology (ext 2436)

VALPROIC ACID, Code: VAL

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: **Therapeutic Drug Monitoring**
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

VANCOMYCIN LEVEL (RANDOM), Code: VAN

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: NA
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: NA
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

VANCOMYCIN LEVEL (PEAK), Code: VANP

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Collect 30 min after end of IV infusion or 60 min after IM injection and send to the Lab immediately
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

VANCOMYCIN LEVEL (TROUGH), Code: VANT

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: Collect specimen immediately before the next dose and send to the Lab immediately
 Specimen Stability: 2 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Particle-enhanced turbidimetric inhibition immunoassay (PETINIA)
 Useful for: Therapeutic Drug Monitoring
 Reference range: See laboratory report
 Performed: Daily
 TAT: 1 day, STAT
 Section: Biochemistry (ext 2407)

Refer to changes highlighted yellow

VITAMIN B12, Code: B12

Specimen: Serum (SST yellow top, 5mL) or plasma (Li-Heparin green top, 4 mL)
 Collection: **NA**
 Specimen Stability: 7 days at 2-8°C
 Rejection Criteria: **NA**
 Method: Chemiluminescent Microparticle Immunoassay (CMIA)
 Useful for: Assessment of vitamin B12 deficiency
 Reference range: See laboratory report
 Performed: Every Tuesday and Friday
 TAT: 3 days
 Section: Immunology (ext 2407)

WIDAL/WEIL FELIX, Code: WWF

Specimen: Serum (SST tube, Yellow top, 2 mL)
 Collection: -
 Specimen Stability: 2°C to 8°C for 48 hours
 Rejection Criteria: Grossly Haemolysed/ Grossly Lipaemic/Contaminated
 Method: Agglutination
 Useful for: Diagnosis of *Salmonella* and Rickettsial infection
 Reference range: Titre < 1:80
 Performed: Office Hours only
 TAT: 1 day
 Section: Microbiology (ext 2436)

WOUND SWAB CULTURE & SENSITIVITY, Code: WSC

Specimen: Wound (Gel Swab/Sterile container/Syringe with needle removed)
 Collection: Indicate source of specimen
 Specimen Stability: 24 hours at room temperature
 Rejection Criteria: Dry swab, samples with formalin
 Method: Conventional culture and Gram Stain
 Useful for: Diagnosis of bacterial infection
 Reference range: -
 Performed: Daily
 TAT: 48-72 hours
 Section: Microbiology (ext 2436)

ANNEX B. LIST OF STAT TESTS

- ◇ STAT (Short-Turn-Around-Time), or Urgent, tests are highest priority tests ordered to manage medical emergencies. The following tests can be performed as STAT (Short Turn Around Time):

BIOCHEMISTRY	Amikacin (AMIK/AMIKP/AMIKT)	1 hour	Neonatal Bilirubin (NBI)	1 hour
	Blood Gas (GAS/VBG)	30 mins	Osmolality (OSS)	2 hours
	Calcium (CA)	1 hour	Phenobarbital (PHNO3)	1 hour
	Carbamazepine (CARB)	1 hour	Phenytoin (PHNY2)	1 hour
	Chloride (CL)	45 mins	Potassium (POT)	45 mins
	Creatinine (CRT)	45 mins	Protein (Urine Random) (UTP)	1 hour
	Creatine Kinase (CK)	1 hour	Sodium (NA)	45 mins
	Digoxin (DIG)	1 hour	Theophylline (THE)	1 hour
	Gentamicin (GEN/GENP/GENT)	1 hour	Urea (UR)	45 mins
	Glucose Fasting (GLF)	1 hour	Uric Acid (SUA)	45 mins
	Lactic Acid (LAC)	1 hour	Valproic Acid (VAL)	1 hour
	Magnesium (MG)	1 hour	Vancomycin (VAN/VANP/VANT)	1 hour
IMMUNOLOGY	N-Terminal Pro B-Type Natriuretic Peptide (NTPROBNP)	1 hour	Methotrexate (MTX)	1 hour
	Beta HCG Total (QUA)	1 hour	Procalcitonin (PCT)	1 hour
	Cortisol (RCO)	1 hour	Progesterone (PRG)	1 hour
	Cyclosporine (CYCLO)	2 hours	Tacrolimus (TACRO)	2 hours
	Estradiol (E2)	1 hour	Troponin I Hs Stat (TRO)	45 mins
HAEMATOLOGY	Activated Partial Thromboplastin Time (APTT)	1 hour	International Normalised Ratio (INR)	1 hour
	D-Dimer (FDP)	2 hours	Malarial Parasites (MS)	1 hour
	Fibrinogen (FIB)	2 hours	Prothrombin Time (PT)	1 hour
	Full Blood Count (FBC)	45 mins	Reticulocyte Count (RET)	45 mins
MICROBIOLOGY	Urine Microscopy (MU)	1 hour	Fluid FEME (FLF)	2 hours
	Faecal Microscopy (FM)	1 hour	CSF FEME (CSFF)	2 hours
VIROLOGY	RSV (< 5 years old)	40 mins	Influenza Antigen	40 mins

Refer to changes highlighted yellow

ANNEX C. TEST CATALOGUE (OUTSOURCE)

For any outsourced test enquiry (including tests not listed below), please call ext. 2127 during office hours only.

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
Acetoacetate (Ketone)	KET	2x 4ml EDTA (whole blood)	30	
Acetylcholinesterase Receptor A	ARA	8ml GEL Yellow Top	30	
Acid Fast Bacilli Smear and Culture	AFBCS	Sputum, CSF, Urine or Body Fluid	60	Provide BruHims No.
Activated Protein C Resistance	APCR	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	59	
Adrenocorticotrophic Hormone (ACTH)	ACT	2 x 4ml EDTA (purple) Frozen Plasma	17	
Mycobacterium Tuberculosis Quantiferon TB Gold	MTBG	Using specific collection kit: (Obtain from Lab)	17	-Specific Request Form -Sample must be collected and send to Lab: Monday until Friday before 09:30AM
Aldolase (Serum)	ALD	8ml GEL Yellow Top	17	
Aldosterone (Serum)	ALDOS	2 x 4ml EDTA (purple) Frozen Plasma	30	
Alpha 1 Antitrypsin	AAT	8ml GEL Yellow Top	30	
Angiotensin Converting Enzyme	ACE	8ml GEL Yellow Top – Frozen Serum	30	
Anti Diuretic Hormone	ADH	2 x 4ml EDTA (purple) Frozen Plasma	30	
Anti Mitochondrial Antibodies	AMA	8ml GEL Yellow Top	30	
Anti Neutrophil Cytoplasmic Antibody, Myeloperoxidase Antibody, Proteinase 3 Antibody	ANCAR	8ml GEL Yellow Top	30	Provide BruHims No.
Anti Nuclear Antibody	ANAR	8ml GEL Yellow Top	30	Provide BruHims No.
Anti Smooth Muscle Antibody	ASM	8ml GEL Yellow Top	17	
Anti Smooth Muscle, Anti Parietal Cell Ab, Anti Mitochondrial Ab	MSKR	8ml GEL Yellow Top	30	Provide BruHims No.

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
Anti Thrombin III	AT3	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	30	
Anti-dsDNA	DNAR	8ml GEL Yellow Top	30	Provide BruHims No.
Anti-Endomysial Ab (IgA)	AEA	8ml GEL Yellow Top	17	
Anti-mullerian hormone	AMH	8ml GEL Yellow Top	10	Preferably Monday until Thursday
Anti-Myelin Oligodendrocyte Glycoprotein Ab	MOG	8ml GEL Yellow Top	10	
Anti-Titin Antibody	T11	8ml GEL Yellow Top Or 1ml CSF	17	
Aquaporin Antibody	NMO	8ml GEL Yellow Top Or 1ml CSF	30	
Babysafe IEM Screening	IEM	Guthrie Card	30	Specific Request Form*
Beta 2 Glycoprotein 1 IgG & IgM	GPITR	8ml GEL Yellow Top	30	Provide BruHims No.
Beta 2 Microglobulin	BMG	8ml GEL Yellow Top	30	
Beta X Laps (C terminal telopeptide)	BXL	8ml GEL Yellow Top – Frozen Serum	30	
Bone Marrow Examination	BMEXAM	Bone marrow	14	Pre-arrange with Lab (ext 2408)
BRCA 1,2, PTEN & p53 Germline(Blood)	BRCAB	2 x 4ml EDTA Blood (Purple)	30	-Specific Request Form -Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Caeruloplasmin	CAE	8ml GEL Yellow Top	17	
Calcitonin	CLT	8ml GEL Yellow Top (frozen serum)	30	
Cardiolipin Antibodies IgG & IgM	CPIR	8ml GEL Yellow Top	30	Provide BruHims No.
CD-19 (B-Cell) Assay	CD19	2 x Fresh 4ml EDTA (Purple Top)	30	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
CD-20 (PAN B) Assay	CD20	2 x Fresh 4ml EDTA (Purple Top)	30	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Chlamydia IgG	CHG	8ml GEL Yellow Top	10	

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
Chromogranin A	CHGA	8ml GEL Yellow Top	17	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Chromosome Analysis (Blood)	CHB	8ml LITHIUM HEPARIN (green)	17	-Specific Request Form -Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Clopidogrel (CYP2C19)	CYP2	2 x 4ml EDTA Blood (Purple)	17	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Copper (24hrs Urine)	UCU	24hr Urine Bottle	17	
C-Peptide	CPP	8ml GEL Yellow Top	17	
Cytomegalovirus IgG	CMG1	8ml GEL Yellow Top	10	Provide BruHims No.
Cytomegalovirus IgM	CMM1	8ml GEL Yellow Top	10	Provide BruHims No.
Cytomegalovirus PCR QUALITATIVE	CMVPCR	Fluid, CSF, 4ml Whole Blood (Purple top), or Dry Swab	17	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Down Syndrome (Triple Test)	NTD	8ml GEL Yellow Top (15w0d – 19w6d)	10	
Drug Confirmation Test	GCMS	20ml Urine Drug	30	
Encephalopathy Autoimmune Evaluation	ENS2	8ml GEL Yellow Top	30	
Epidermal Growth Factor Receptor Mutation (Blood)	EGFR	2x8ml cell streak tube (Obtain from Lab)	30	-Specific Request Form -Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Epstein Barr Virus Ab VCA IgG	EBG	8ml GEL Yellow Top	17	
Epstein Barr Virus Ab VCA IgM	EBM	8ml GEL Yellow Top	17	
Epstein Barr Virus EA IgA (NPC)	EBE	8ml GEL Yellow Top	10	
Epstein Barr Virus PCR	EBVPCR	2 x Fresh 4ml EDTA (Purple Top)	17	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Epstein Barr Virus VCA IgA (NPC Marker)	EBA	8ml GEL Yellow Top	10	

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
Extractable Nuclear Antigen	ENA	8ml GEL Yellow Top	17	
Factor 8	F8	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	30	
Factor 9	F9	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	30	
First Trimester Screen	FTS	8ml GEL Yellow Top	10	Specific Request Form
Flavivirus RT PCR	FLAVI	2 x 4ml EDTA (purple) Frozen Plasma; OR 1 X CSF Liquid	59	
Foundation One Liquid	FM1L	2 x 8ml Streck Tube	30	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Free PSA & Total PSA	FCA	8ml GEL Yellow Top	17	
Gangliosides Profile (GM1, GM2, GM3, GD1a, GD1b, GT1b, GQ1b)	GANG	2ml CSF	17	
Gangliosides Profile (GM1, GM2, GM3, GD1a, GD1b, GT1b, GQ1b)	GM2B	8ml GEL Yellow Top	30	
Glutamic Acid Decarboxylase Ab	GAD	8ml GEL Yellow Top or 8ml Red Top	30	
Growth Hormone	GH	8ml GEL Yellow Top – Frozen Serum	17	
Haemoglobin Electrophoresis	HBE	8ml GEL Yellow Top, AND 2 x 4ml EDTA(Purple) Whole Blood	10	
Haptoglobin	HAP	8ml GEL Yellow Top	17	
Helicobacter Breath Test	HEB	Urea Beath Test Kit (Obtain from Lab)	10	
Helicobacter Pylori IgG	HEL	8ml GEL Yellow Top	10	
Hepatitis B Core IgM	HBC	8ml GEL Yellow Top	7	Provide BruHims No.
Hepatitis B Envelope Antibody	HBV	8ml GEL Yellow Top	7	Provide BruHims No.
Hepatitis B Envelope Antigen	HBX	8ml GEL Yellow Top	7	Provide BruHims No.
Hepatitis C Antibody Confirmation	HCVC	8ml GEL Yellow Top	7	Provide BruHims No.
Hepatitis C RNA (PCR)	HCQ1	2 x 4ml EDTA (purple) Plasma	14-21	Provide BruHims No.

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
Hepatitis E Antibodies IgG	HEG	8ml GEL Yellow Top	17	
Hepatitis E Antibodies IgM	HEM	8ml GEL Yellow Top	17	
Herpes Simplex I IgG	HG1	8ml GEL Yellow Top	17	
Herpes Simplex II IgG	HG2	8ml GEL Yellow Top	17	
Herpes Simplex Type 1 and 2 (IgM)	HM12	8ml GEL Yellow Top	17	
Histology	HIR	Container with 10% buffered formalin	14	Provide BruHims No.
Histology (More than 7 bottles)	HIR2	Container with 10% buffered formalin	14	Provide BruHims No.
HIV-1/2 Antibodies Confirmation	HIC	8ml GEL Yellow Top	7	Provide BruHims No.
HRR Blood	HRRB	2 x Fresh 4ml EDTA (Purple Top)	30	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Human Papillomavirus PCR	HPVPCR	Using specific kit (get the kit from Laboratory)	17	
Immunoglobulin A (IgA)	IGA	8ml GEL Yellow Top	17	
Immunoglobulin E (IgE)	IGE	8ml GEL Yellow Top	17	
Immunoglobulin IgG, IgM & IgA	EPP	8ml GEL Yellow Top	10	
Immunoglobulin M (IgM)	IGM	8ml GEL Yellow Top	30	
Immunoglobulin Subclass IgG4	IGGS4	8ml GEL Yellow Top	30	
Insulin Like Growth Factor (IGF1-Stomatomedic C)	IGF	8ml GEL Yellow Top	30	
Intrinsic Factor Antibody	IFA	8ml GEL Yellow Top	17	
JAK 2 Mutation	JAK	2 x Fresh 4ml EDTA (Purple Top)	30	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
John Cunningham Virus (JCV) PCR	JCVP	2 x Fresh 4ml EDTA (Purple Top) Whole Blood OR 2ml CSF	30	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
KRAS & NRAS Mutation (Blood)	KNRAS	2x8ml cell streak tube (Obtain from Lab)	30	-Specific Request Form -Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Lamotrigine	LAMO	8ml GEL Yellow Top	30	
Legionella Antigen	LEGU	20ml Random Urine	17	
Leptospirosis Antibody (IgM)	LEPM	8ml GEL Yellow Top	17	
Leptospirosis IgG	LEPG	8ml GEL Yellow Top	10	
Levetiracetam (Keppra Level)	KEP	8ml GEL Yellow Top	17	
Lipase	LSE	8ml GEL Yellow Top	17	
Lithium	LIT	8ml GEL Yellow Top	10	
Lupus Anti-Coagulant Factor	LAS	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	30	
Measles IgG Antibody	MEG1	8ml GEL Yellow Top	10	Provide BruHims No.
Measles IgM Antibody	MEM	8ml GEL Yellow Top	10	Provide BruHims No.
Melioidosis Antibody IgG	MELG	8ml GEL Yellow Top	30	
Melioidosis Antibody IgM	MELM	8ml GEL Yellow Top	30	
Metanephrine in 24 Hours Urine	MNP	24hr Urine bottle with 30ml 6N HCL	30	
Mumps IgG	MUMG1	8ml GEL Yellow Top	10	Provide BruHims No.
Mumps IgM	MUMM	8ml GEL Yellow Top	10	Provide BruHims No.
Mycoplasma Pneumonia IGM	MYCO	8ml GEL Yellow Top	10	
Mycoplasma Pneumonia Total AB	MYCOAB	8ml GEL Yellow Top	10	
Myoglobin (Urine)	UMYO	20ml random urine	17	
Myositis Associated Antibody	MA	8ml GEL Yellow Top	17	
Neuron Specific Enolase	NSE	8ml GEL Yellow Top	30	
Neuronal Antibodies Extended Panel (12 Abs)	NAEP	8ml GEL Yellow Top Or CSF	30	

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
NEURONAL Antigen Profile Immunoblot	NEUROBLOT	8ml GEL Yellow Top	30	Provide BruHims No.
NICE BASIC (Non-invasive Prenatal Test)	NICEB	Using specific collection kit: (get the kit from laboratory)	17	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
NICE EXTENDED (Non-invasive Prenatal Test)	NICEX	Using specific collection kit: (get the kit from laboratory)	17	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
NMDA Ab	NMDA	8ml GEL Yellow Top Or 1ml CSF	30	
Paraneoplastic Antibodies	PAR	8ml GEL Yellow Top	10	
Personal Cancer Diagnostic Invtiae	PCDI	2 x 4ml EDTA Blood (Purple)	30	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Phospholipid Antibody Screen	CPILAS	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma, AND 1 x 8ml GEL Yellow Top	17	
Pik3Ca Liquid Biopsy	LBPIK3CA	2 x Streck tube (Obtain from lab)	17	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Protein C	PTC	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	17	
Protein Electrophoresis	PES	8ml GEL Yellow Top	30	
Protein S	PTS	2 x 2.7ml Blue Top (Sodium Citrate) Frozen Plasma	17	
Prothrombin Mutation	PTG	2 x 4ml EDTA Blood (Purple)	30	- Specific Request Form - Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Renin	REN	2 x 4ml EDTA (purple) -Frozen Plasma	17	

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
Salicylates	SAL	8ml GEL Yellow Top	30	
SHBG (Sex Hormone Binding Globulin)	SHB	8ml GEL Yellow Top	17	
Stone (Calculi) Analysis	STO	Stone in sterile container	14	
SurePath PAP Smear	TPP	Using specific kit	14	
Testosterone (Free)	TESF	8ml GEL Yellow Top	17	Sample must be collected and send to Lab – Monday, Tuesday or Thursday before 09:30AM
Thalassaemia Studies	THS	8ml GEL Yellow Top, AND 2 x 4ml EDTA (purple)	17	
Thiamine (Vitamin B1)	VB1	2 x 4ml EDTA Whole Blood (Purple top), wrap with Aluminium foil	30	
Total T3	T3	8ml GEL Yellow Top	10	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Total T4	T4	8ml GEL Yellow Top	10	Sample must be collected and send to Lab: Monday, Tuesday or Thursday before 09:30AM
Toxoplasma IgG	TXG	8ml GEL Yellow Top	17	
Toxoplasma IgM	TXM	8ml GEL Yellow Top	17	
Transglutaminase Antibody (IgA & IgG)	TRANAB	8ml GEL Yellow Top	30	
Tuberculosis PCR	TBPCR	Sputum, CSF, Urine or Body Fluid	30	Provide BruHims No.
Urine Catecholamines 24 Hrs	UCAT	24hr Urine Bottle with acid (30ml 6N HCL added)	17	
Urine Cytology	UCYT	20ml random urine	17	Provide BruHims No.
Urine Porphyrin Profile	UPSG	20ml random urine (wrap in foil and freeze)	17	
Urine Protein Electrophoresis	UPES	20ml random urine	17	
Urine Vanillylmandelic Acid 24	VMA	24hr Urine Bottle in 20ml 6N HCL	17	
Varicella Zoster IgG	VZG	8ml GEL Yellow Top	7	Provide BruHims No.
Varicella Zoster IgM	VZM	8ml GEL Yellow Top	7	Provide BruHims No.

JPMC Laboratory Handbook

Refer to changes highlighted yellow

Test Name	Test Code	Specimen Requirement	TAT (working days)	Remarks
VDRL for CSF	VDRLCSF	2ml CSF Fluid	10	
Ziehl-Neelsen stain	ZN	Any (Sterile container)	7	

CALL OUTSOURCE SECTION (ext 2127) FOR SPECIFIC REQUEST FORM AND SPECIMEN REQUIREMENTS

ANNEX D. TEST PACKAGES

JPMC PACKAGES

**Packages are subject to change. Please call OPD for availability of test packages.*

JPMC PREMIER HEALTH SCREENING PACKAGE, Code: JHSP
 EMPLOYEE HEALTH PACKAGE, Code: LEHP
 JPMC EXCLUSIVE MALE HEALTH SCREENING PACKAGE, Code: JHSPM
 JPMC EXCLUSIVE FEMALE HEALTH SCREENING PACKAGE, Code: JHSPEF
 GENERAL TEST 1 PACKAGE, Code: ZGT1
 GENERAL TEST 5 PACKAGE, Code: ZGT5
 GENERAL TEST 7 PACKAGE, Code: ZGT7
 GENERAL TEST FOR PAEDIATRICS, Code: ZGTP
 ANTENATAL PACKAGE, Code: ZANJ [OBYN CLINIC]

JPMC PREMIER HEALTH SCREENING PACKAGE, Code: JHSP

Full Blood Examination (including ESR)	1 X PURPLE TUBE	WHOLE BLOOD
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
Hepatitis B Screen (HBsAg + HBsAb)		
Hepatitis C Total Antibody		
Fasting or Random Glucose	1 x GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE

EMPLOYEE HEALTH PACKAGE, Code: LEHP

Malaria Screen	1 X PURPLE TUBE	WHOLE BLOOD
Hepatitis B Screen (HBsAg & HBsAb)	1 X YELLOW TUBE	SERUM
HIV 1/2 Ag/Ab		

JPMC Laboratory Handbook

Refer to changes highlighted yellow

JPMC EXCLUSIVE MALE HEALTH SCREENING PACKAGE, Code: JHSPeM

Full Blood Examination (including ESR)	1 X PURPLE TUBE	WHOLE BLOOD
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
Alpha-Feto Protein		
Hepatitis B Screen (HBsAg + HBsAb)		
Hepatitis C Total Antibody		
Cancer Marker 19.9		
Carcinoembryonic Antigen		
Total Prostate Specific Antigen		
Epstein Barr Virus EA IgA (NPC Marker) Send out		
Fasting Glucose	1 X GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE
Faecal Occult Blood	STERILE CONTAINER	STOOL

JPMC EXCLUSIVE FEMALE HEALTH SCREENING PACKAGE, Code: JHSPEF

Full Blood Examination (including ESR)	1 X PURPLE TUBE	WHOLE BLOOD
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
Alpha-Feto Protein		
Hepatitis B Screen (HBsAg + HBsAb)		
Hepatitis C Total Antibody		
Cancer Marker 19.9		
Cancer Marker 125		
Cancer Marker 15.3		
Carcinoembryonic Antigen		
Fasting Glucose	1 X GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE
Faecal Occult Blood	STERILE CONTAINER	STOOL
SurePath Pap Test Send out		

JPMC Laboratory Handbook

Refer to changes highlighted yellow

GENERAL TEST 1 PACKAGE, Code: ZGT1

Full Blood Examination (including ESR)	1 X PURPLE TUBE	WHOLE BLOOD
Renal Function Test	1 X YELLOW TUBE	SERUM
Liver Function Test		
Lipid Profile		
Fasting Glucose	1 x GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE

GENERAL TEST 5 PACKAGE, Code: ZGT5

Full Blood Count (including ESR)	2 X PURPLE TUBES	WHOLE BLOOD
Blood Group		
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
Hepatitis B Screen (HBsAg + HBsAb)		
Rheumatoid Factor		
Hepatitis A Antibodies (IgM + IgG)		
Free Thyroxine (FT4)		
Alpha-Feto Protein		
Carcinoembryonic Antigen		
Fasting Glucose		
Urine FEME	STERILE CONTAINER	RANDOM URINE

JPMC Laboratory Handbook

Refer to changes highlighted yellow

GENERAL TEST 7 PACKAGE, Code: ZGT7

Full Blood Count (including ESR)	2 X PURPLE TUBES	WHOLE BLOOD
Blood Group		
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
Hepatitis B Screen (HBsAg + HBsAb)		
Rheumatoid Factor		
Hepatitis A Antibodies (IgM + IgG)		
Free Thyroxine (FT4)		
Alpha-Feto Protein		
Total Prostate Specific Antigen		
Fasting Glucose	1 X GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE

GENERAL TEST FOR PAEDIATRICS, Code: ZGTP

Full Blood Count (with Blood Film)	1 X PURPLE TUBE (or 1 X PURPLE MICROTAINER)	WHOLE BLOOD
Renal Function Test	1 X YELLOW TUBE (or 2 X GREEN MICROTAINERS)	SERUM/ PLASMA
Liver Function Test		
Calcium		
Phosphate		
C-Reactive Protein		
Fasting Glucose		

ANTENATAL PACKAGE, Code: ZANJ [OBGYN CLINIC]

Full Blood Count	2 X PURPLE TUBES	WHOLE BLOOD
Blood Group		
Antibody Screen (Indirect Coombs Test)		
Hepatitis B Screen (HBsAg + HBsAb)	1 x YELLOW TUBE	SERUM
Rubella IgG Antibodies		
Random Glucose	1 X GREY TUBE	PLASMA

JPMC Laboratory Handbook

Refer to changes highlighted yellow

BNSRC PACKAGES

**Only available for BNSRC Clinic and Wards*

BNSRC PROFILE TEST 1, Code: NSRC1

BNSRC PROFILE TEST 2, Code: NSRC2

BNSRC PROFILE TEST 3, Code: NSRC3

BNSRC PROFILE TEST 1, Code: NSRC1

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Urea	1 X YELLOW TUBE	SERUM
Electrolytes		
Creatinine		
Calcium		
Total Protein		
Gamma-Glutamyl Transferase		
Alanine Aminotransferase		
Aspartate Aminotransferase		
Creatine Kinase		
C-Reactive Protein		
Lactate Dehydrogenase		
Thyroid Stimulating Hormone		
Troponin I		
Coagulation Profile	1 X LIGHT BLUE TUBE	PLASMA
Glucose Random	1 X GREY	PLASMA

BNSRC PROFILE TEST 2, Code: NSRC2

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes	1 X YELLOW TUBE	SERUM
C-Reactive Protein		
Creatine Kinase		
Creatinine		
Gamma-Glutamyl Transferase		
Alanine Aminotransferase		
Aspartate Aminotransferase		

JPMC Laboratory Handbook

Refer to changes highlighted yellow

BNSRC PROFILE TEST 3, Code: NSRC3

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Urea	1 X YELLOW TUBE	SERUM
Electrolytes		
C-Reactive Protein		
Creatinine		
Creatinine Kinase		
Gamma-Glutamyl Transferase		
Alanine Amino Transferase		
Aspartate Aminotransferase		
Calcium		
Albumin		
Coagulation Profile	1 X LIGHT BLUE TUBE	PLASMA

TBCC PACKAGES

**Only available for TBCC Clinic and Wards*

BREAST CANCER PACKAGE, Code: BRCP
 COLORECTAL & GIT CANCER PACKAGE, Code: COCP
 GERM CELL TUMOUR PACKAGE, Code: GCCP
 GYNAECOLOGICAL CANCER PACKAGE, Code: GYCP
 HEAD & NECK CANCER PACKAGE, Code: HNCP
 LIVER CANCER PACKAGE, Code: LICP
 LUNG CANCER PACKAGE, Code: LUCP
 LYMPHOMA CANCER PACKAGE, Code: LYMCP
 NCC PACKAGE, Code: NCCP
 NPC PRE-TREATMENT PACKAGE, Code: NPC1
 NPC DURING CHEMOTHERAPY PACKAGE, Code: NPC2
 PANCREAS AND BILIARY CANCER PACKAGE, Code: PACP
 PROSTATE CANCER PACKAGE, Code: PRCP
 THYROID PRE-TREATMENT, Code: TPTP1
 THYROID DURING CHEMOTHERAPY, Code: TPTP2

BREAST CANCER PACKAGE, Code: BRCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Cancer Marker CA 15-3		
Glucose Random	1 X GREY TUBE	PLASMA

COLORECTAL & GIT CANCER PACKAGE, Code: COCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Carcinoembryonic Antigen		
Glucose Random	1 X GREY TUBE	PLASMA

JPMC Laboratory Handbook

Refer to changes highlighted yellow

GERM CELL TUMOUR PACKAGE, Code: GCCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Alpha-Feto Protein		
Beta Human Chronic Gonadotrophin		
Lactate Dehydrogenase		
Glucose Random	1 X GREY TUBE	PLASMA

GYNAECOLOGICAL CANCER PACKAGE, Code: GYCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Cancer Marker CA 12-5		
Glucose Random	1 X GREY TUBE	PLASMA

HEAD & NECK CANCER PACKAGE, Code: HNPC

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Glucose Random	1 X GREY TUBE	PLASMA

LIVER CANCER PACKAGE, Code: LICP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Alpha-Feto Protein		
International Normalised Ration (INR)	1 X LIGHT BLUE TUBE	PLASMA
Glucose Random	1 X GREY TUBE	PLASMA

JPMC Laboratory Handbook

Refer to changes highlighted yellow

LUNG CANCER PACKAGE, Code: LUCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Carcinoembryonic Antigen		
Cancer Marker CA 19-9		
Glucose Random	1 X GREY TUBE	PLASMA

LYMPHOMA CANCER PACKAGE, Code: LYMCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Erythrocyte Sedimentation Rate	1 X ESR TUBE	
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM/PLASMA
Liver Function Test		
Lactate Dehydrogenase		
Glucose Random	1 X GREY TUBE	PLASMA

NCC PACKAGE, Code: NCCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Glucose Random	1 X GREY TUBE	PLASMA

NPC PRE-TREATMENT PACKAGE, Code: NPC1

Note: ZNPC1 specimen must be taken from Monday, Tuesday and Thursday (before 09:30 AM) only

Full Blood Count	2 X PURPLE TUBES	WHOLE BLOOD
Epstein Barr Virus PCR (Send out)		
Electrolytes, Creatinine & Urea	2 X YELLOW TUBES	SERUM
Liver Function Test		
Epstein Barr Virus VCA IgA (Send out)		
Glucose Random	1 X GREY TUBE	PLASMA

JPMC Laboratory Handbook

Refer to changes highlighted yellow

NPC DURING CHEMOTHERAPY PACKAGE, Code: NPC2

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Glucose Random	1 X GREY TUBE	PLASMA

PANCREAS AND BILIARY CANCER PACKAGE, Code: PACP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Cancer Marker CA 19-9		
Glucose Random	1 X GREY TUBE	PLASMA

PROSTATE CANCER PACKAGE, Code: PRCP

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Total Prostate Specific Antigen		
Glucose Random	1 X GREY TUBE	PLASMA

THYROID PRE-TREATMENT, Code: TPTP1

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	2 X YELLOW TUBES	SERUM
Liver Function Test		
Thyroid Function Test		
Thyroglobulin (Send out)		
Glucose Random	1 X GREY TUBE	PLASMA

THYROID DURING CHEMOTHERAPY, Code: TPTP2

Full Blood Count	1 X PURPLE TUBE	WHOLE BLOOD
Electrolytes, Creatinine & Urea	1 X YELLOW TUBE	SERUM
Liver Function Test		
Glucose Random	1 X GREY TUBE	PLASMA

JPMC Laboratory Handbook

Refer to changes highlighted yellow

GJPMC PACKAGES

**Only available for GJPMC Clinic and Wards*

GENERAL TEST 1 PACKAGE, Code: GT1GJ
 GENERAL TEST 5 PACKAGE, Code: GT5GJ
 GENERAL TEST 7 PACKAGE, Code: GT7GJ
 OPEN HEART SURGERY SCREENING, Code: OHS

GENERAL TEST 1 (GJ) PACKAGE, Code: GT1GJ

Full blood count (including ESR)	2 X PURPLE TUBES	WHOLE BLOOD
Glycated Haemoglobin (HbA1c)		
Renal function test	1 X YELLOW TUBE	SERUM
Liver function test		
Lipid Profile		
Glucose fasting	1 X GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE

GENERAL TEST 5 (GJ) PACKAGE, Code: GT5GJ

Full blood count (including ESR)	3 X PURPLE TUBES	WHOLE BLOOD
Blood Group		
Glycated Haemoglobin (HbA1c)		
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
Hepatitis B Screen (HBsAg + HBsAb)		
Rheumatoid Factor		
Hepatitis A Antibodies (IgM + IgG)		
Thyroid Function Test		
Alpha Feto Protein		
Carcinoembryonic Antigen		
Cancer Marker 19.9		
Cancer Marker 125		
Cancer Marker 15.3		
RPR Screen		
Glucose Fasting	1 X GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE

JPMC Laboratory Handbook

Refer to changes highlighted yellow

GENERAL TEST 7 (GJ) PACKAGE, Code: GT7GJ

Full Blood Count (Including ESR)	3 X PURPLE TUBES	WHOLE BLOOD
Blood Group		
Glycated Hemoglobin (HbA1c)		
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver function test		
Lipid Profile		
Hepatitis B Screen (HBsAg + HBsAb)		
Rheumatoid Factor		
Hepatitis A Antibodies (IgM + IgG)		
Thyroid Function Test		
Alpha Feto Protein		
Carcinoembryonic Antigen		
Cancer Marker 19.9		
Total Prostate Specific Antigen		
RPR Screen		
Glucose Fasting	1 X GREY TUBE	PLASMA
Urine FEME	STERILE CONTAINER	RANDOM URINE

OPEN HEART SURGERY SCREENING, Code: OHS

Full Blood Examination (Including ESR)	3 X PURPLE TUBES	WHOLE BLOOD
Blood Group And Hold		
Glycated Hemoglobin (HbA1c)		
Renal Function Test	2 X YELLOW TUBES	SERUM
Liver Function Test		
Lipid Profile		
C-Reactive Protein		
RPR Screen		
Hepatitis B Screen (HBsAg + HBsAb)		
HIV 1 & 2 Ag/Ab		
Hepatitis A Antibodies (IgM + IgG)		
Hepatitis C Total Antibody		
Thyroid Function Test		
Glucose Fasting	1 X GREY TUBE	PLASMA
Coagulation Profile	1 X LIGHT BLUE TUBE	PLASMA
Urine Culture And Sensitivity	STERILE CONTAINER	RANDOM URINE
MRSA Swab For Culture And Sensitivity	3 X SWABS	GROIN, AXILLA & NASAL
Sterility Test		
Sputum Culture And Sensitivity	STERILE CONTAINER	SPUTUM